Laparoscopic ABC® Probes™ - Single Patient Use

- REF 130342 - 10mm Footswitching (28cm)
- REF 160655 - 5mm Footswitching (28cm)
- REF 160656 - 5mm Handswitching (28cm)
- REF 160636 - 5mm Handswitching (36cm)
- REF 160644 - 5mm Handswitching (44cm)

Sterility guaranteed unless package has been opened, damaged or compromised.

Description:

The Laparoscopic Argon Beam Coagulation (ABC®) Extended Probe is a single patient use product, provided sterile. This argon gas enhanced monopolar electrosurgical device is intended to be used through a laparoscopic trocar sleeve.

Indications:

Use the Laparoscopic ABC® Extended Probe only in situations where you would normally use monopolar energy.

For use with CONMED® Electrosurgical Generators with ABC® and ABC® Modules.

Maximum generator voltage 6.5kV peak @ 3% duty cycle.

Contraindications

- Do not use this device in situations where conventional monopolar electrosurgery is inappropriate or unsafe.
- Do not use this argon gas enhanced device in procedures where the risk of gas embolism exists (such as inside the uterus).
- Do not direct the nozzle into open vessels or bury the nozzle into tissue. Intra-abdominal pressure should not exceed 18 millimeters of Hg. Discontinue use of this device if intra-abdominal pressure exceeds this value. Failure to follow these warnings and instructions (i.e. maximum flow rate of 4 liters per minute of argon and 18mm Hg intra-abdominal pressure) could increase risk of embolism.
- Do not use in the presence of flammable anesthetics, disinfecting agents, oxygen-rich environments or other combustible materials.

Warnings:

- This device should be used only in conjunction with an electronic insufflator and venting trocar to control the intra-abdominal pressure.
- This is a monopolar, active instrument. Before use, a dispersive electrode (i.e. grounding pad) must be applied to the patient as recommended by the manufacturer.
- Do not directly secure the handpiece cord to surgical drapes with metal instruments. Activation of the handpiece in contact with metal instruments may cause burns at the tissue/instrument interface.
- Isolate the handpiece in an accessory holster when not in use. This will prevent inadvertent contact with fluids, the patient, or surgical personnel.
- Question unusual requests for increased power output settings. Such requests should initiate a check of integrity/condition of all cords, connections, and electrodes including the dispersive electrode (ground pad) before a power increase is made.
- Care should be taken during surgery to prevent inadvertent contact of the handpiece/probe tip with tissue.
- Do not activate the electrosurgical unit, if the tip of the instrument is not in a position to deliver energy (fulguration) to target tissue. Doing so may cause inadvertent burns to patient.
- **MODIFICATION OF THE DEVICE BY THE USER:** CONMED® warns explicitly against modifying the device, e.g., cutting off the tip of the malleable tube. Any modification voids warranty and exempts CONMED® from all liability.
- Electrosurgery is potentially hazardous for patients with pacemakers, AICDs, neurostimulators or other active implants. High frequency energy may damage or adversely affect active implants. Ventricular fibrillation may occur.

Precautions:

- Refer to CONMED® ABC® Generator operator manual and CONMED® ABC® Handpiece instructions prior to use.
- **SINGLE PATIENT USE. DO NOT RESTERILIZE.** The ability to effectively clean and re-sterilize this device has not been established and subsequent re-use may adversely affect the performance, safety and/or sterility of the device. Improper operation of this device may compromise device safety.
- Use of the handpiece without gas flow may cause severe tip damage. Make certain the ABC® connector is fully seated into the gas supply receptacle and make sure the gas supply hose does not become kinked or crimped. (Always test for gas flow prior to surgical use by activating the handpiece and directing the active tip toward pooled fluid. Appropriate gas flow will result in fluid dispersion prior to establishment of Argon Beam.)
- Before use, inspect the cord insulation and handpiece integrity and condition. If damaged, chipped, cut, or nicked, do not use the handpiece/probe.
- Do not immerse, or pour fluids over, the handpiece.
- Do not pull or stretch the cord.
- Keep the handpiece and cord away from sharp objects that may cause damage to insulation.
- Do not test the generator by sparking the active electrode to the dispersive electrode (ground pad) or other objects.
- Argon Beam Coagulation should be used only by surgeons experienced in, and aware of, Argon Beam techniques and safety.
- A complete understanding of the principles of monopolar electrosurgical technique is necessary to avoid shocks and burns and the potential for gas embolism to the patient.

System Set-Up Procedure:

1. Use aseptic technique when removing the Handpiece from the package.

2. Remove tip protector from end of probe prior to use.
 3. Inspect the probe for any damage. Do not use if you suspect any damage is present. Notify the manufacturer immediately.
 4. Attach the Laparoscopic ABC® Probe to the ABC® connector and (if using a hand-control Probe being used) the Beam/monopolar hand-control connector on the front panel of the CONMED Electrosurgical Generator with Argon Beam coagulation (integrated units) or ABC® module.
 5. Turn on the main power switch.
 6. Select "Endo" gas flow.
- Warning:**
Gas Flow Rate should not exceed 4 LPM during laparoscopic procedures.
7. Select desired ABC® power level.
 8. Set Gas Flow Rate to 4 liters per minute (LPM).

Recommended for use with a disposable trocar, in order to minimize potential probe damage.

Handpiece/Probe Activation:

Foot Control:

To activate Argon Beam Coagulation, depress ABC® footswitch.

Hand Control:

To activate Argon Beam Coagulation, depress the blue button.

Note:

Handswitching Probes when used with an ABC® footswitch (REF 130146), can also be foot-switch activated.

Handpiece/Probe Operation:

1. Insert ABC® probe through the appropriate trocar by manually opening the valve and inserting the probe.
2. Actively vent in order to prevent high insufflation pressures in the abdominal cavity during laparoscopic use of this device, completely open the trocar stopcock on the trocar containing the ABC® probe prior to and during the activation of the ABC® probe.
3. Depress ABC® button on handpiece or depress ABC® footswitch. Argon gas will begin to flow.
4. Initiation (firing) of the Argon Beam Coagulator will occur as the Probe is advanced to within 1 to 3 centimeters of the tissue.
5. Continuously monitor the abdominal pressure and the patient during use of this device. Do not rely exclusively on the trocar stopcock to maintain the intra-abdominal pressure below 18mm Hg. Every effort should be made to not exceed the 18mm Hg maximum pressure. If, except for very short transient events, the pressure exceeds this value, discontinue use of the Laparoscopic ABC® Probe and allow excess gas to escape. Once pressure has dropped below 18mm Hg you may continue use of the probe.
6. Move the activated Probe within 1-3cm slowly back and forth across the tissue surface in a gentle motion.
7. When adequate hemostasis or tissue effect is achieved, deactivate the beam by releasing the switch, then remove the Laparoscopic ABC® Probe and close the trocar stopcock.

Important:

Arc initiation is dependent on several variables such as ESU type, energy level used, handpiece angle to the target area, and tissue type. It is recommended that the handpiece be held at a 45 degree angle to the tissue for optimum performance and visibility. To initiate the arc, it is recommended that the handpiece tip be brought to a distance within 1cm to the patient tissue, once initiated the handpiece tip may be retracted for a more comfortable position and arc length.

Do not touch the nozzle directly to tissue or otherwise bury the nozzle tip in tissue or fluids which obscure visualization of the handpiece/probe tip.

After Use:

1. When finished using the handpiece/probe, disconnect it by unplugging the handpiece connectors from the generator and argon gas unit.
2. After disconnection, this product may be a potential biohazard. Handle and disperse of the handpiece as required by hospital policy and applicable law. Do not re-sterilize.

CAUTION: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

525 French Road • Utica, NY 13502-5994 USA
 (315) 797-8375 • FAX (315) 797-0321
 Customer Service (USA) 1-800-448-6506
 e-mail: info@conmed.com
 http://www.conmed.com

European Authorized Representative
 MDSS GmbH
 Schiffgraben 41
 D - 30175 Hannover
 Germany
 For International orders or inquiries,
 please contact
 CONMED International Sales
 (315) 797-8375 • FAX (315) 735-6235

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© CONMED Corporation, Imprimé aux États-Unis d'Amérique