



6875 CHAIR SCALE OPERATING INSTRUCTIONS



WARNING - This product is not intended for use in transporting patients! It is to be used **ONLY** for weighing a patient. Failure to observe this warning may result in serious injury to the patient and/or the scale operator.

1. Install six (6) "AA" size batteries or if ordered with scale, plug the small connector end of the AC power adapter cord into the power jack located in the bottom right corner at the back of the indicator, then plug the AC power adapter into an appropriate power receptacle. **NOTE!** DO NOT connect an AC power adapter to the indicator if Alkaline batteries are installed.
2. Install the indicator on the scale bracket and then connect the load cell cable.

NOTE! For more information on installing the batteries (or the AC power adapter), mounting the indicator on the scale bracket and connecting the load cell cable, refer to the 750 Owner's Manual, 8555-M260-O1.

3. Place the chair scale on any hard, level, flat surface or low-cut carpet.
4. With no weight on scale, press [**ON/OFF**] key. After a second or so, the display should indicate a zero weight with the ZERO annunciator turned on. If display does not indicate zero, press [**ZERO**] key to zero weight display.
5. Set patient on scale and have them remain motionless for several seconds until scale display shows their weight.

IMPORTANT!

**Lower footrest after patient is seated.
Raise footrest before patient stands up.
NEVER STAND UPRIGHT ON FOOTREST!**

6. To re-zero scale, remove patient from scale and press [**ZERO**] key.



0065-C238-1A Rev H
02/10

CARDINAL SCALE MFG. CO.
PO BOX 151 • WEBB CITY, MO 64870
PH (417) 673-4631 • FAX (417) 673-5001
www.detectoscale.com

Printed in USA

Technical Support: Ph: 866-254-8261 • techsupport@cardet.com