Handy-Ox

DIGITAL COLOR PULSE OXIMETER

User Manual

OXÍMETRO DE PULSO DIGITAL

Instrucciones para su uso

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driveTM

INSTRUCTIONS FOR USE

This manual is written in compliance with the council directive MDD93/42/EEC for medical devices and harmonized standards. In case of modifications and software upgrades, the information contained in this document is subject to change without notice.

The manual describes, in accordance with the Digital Handy-Ox Pulse Oximeter's features and requirements, main structure, functions, specifications, correct methods for transportation, installation, usage, operation, repair, maintenance and storage, as well as the safety procedures to protect both the user and equipment. Refer to the respective chapters for details.

Please read the User Manual carefully before using this product. The User Manual describes the operating procedures which should be followed strictly. Failure to follow the User Manual may cause measurement inaccuracies, equipment damage and/or human injury. The manufacturer is NOT responsible for the safety, reliability and performance issues and any monitoring abnormality, human injury and equipment damage due to users' negligence of the operation instructions. The manufacturer's warranty service does not cover such faults.

WARNING:

- An uncomfortable or painful feeling may occur if using the device continuously, especially for the microcirculation barrier patients. It is recommended that the sensor should not be applied to the same finger for over 2 hours.
- For special patients, there should be prudent inspection of the skin when placing the finger probe onto the patient. This device should not be clipped onto an edema or onto tender tissue.
- The light (the infrared is invisible) emitted from the device is harmful to the eyes, so the user and the technician should not stare at the light.
- Patient should not use enamel or other makeup.
- Patient's fingernail should not be too long.
- Please refer to the user instruction manual for clinical use restrictions and cautions.
- This device is not intended for treatment.

Caution: Federal (US) Law restricts this device to sale by or on the order of a physician.

CONTENT

1	Safety	6
	1.1 Instructions for Safe Operations.....	6
	1.2 Warning	6
	1.3 Attention	7
2	Overview	9
	2.1 Features	9
	2.2 Environment Requirements.....	10
3	Principle	10
4	Technical Specifications	11
	4.1 Main Performance.....	11
	4.2 Main Parameters.....	12
5	Installation	14
	5.1 View of the Front Panel.....	14
	5.2 Underside View and Left View.....	14
	5.3 Rear View	15
	5.4 Accessories.....	15
6	Operating Guide	16
	6.1 Application Method.....	16
	6.2 Attention for Operation	26
	6.3 Clinical Restrictions.....	27

7	Maintain, Transportation and Storage	28
	7.1 Cleaning and Disinfecting	28
	7.2 Maintain	28
	7.3 Transportation and Storage.....	28
8	Troubleshooting	29
9	Key of Symbols	30
10	Function Specification	31
	Appendix 1	33
	Appendix 2	37

1 SAFETY

1.1 Instructions for Safe Operations

- Check the main unit and all accessories periodically to make sure that there is no visible damage that may affect patient's safety and monitoring performance. It is recommended that the device is inspected at least once a week. When there is obvious damage, stop using the device.
- Necessary maintenance must be performed by qualified service technicians ONLY. Users are not permitted to maintain the pulse oximeter themselves.
- The oximeter cannot be used together with accessories not specified in User Manual. Only the accessory that is recommended by Drive Medical can be used with this device.
- This product is calibrated before leaving factory.

1.2 Warning

- Explosive hazard—DO NOT use the oximeter in the environment with flammable material such as anesthetic.
- DO NOT use the oximeter while the patient is being scanned by MRI or CT.
- The person who is allergic to rubber cannot use this device.
- The disposal the pulse oximeter and its accessories and packing (including battery, plastic bags, foams and paper boxes) should follow the local laws and regulations.

- Please check the packing before use to make sure the device and accessories are in accordance with the packing list.
- Please only use accessories and probes approved by Drive Medical, or else it may damage the device.
- The device can only be matched with the compatible probe.
- Please don't measure this device with functional tester for the device's related information.

1.3 Attention

- Keep the oximeter away from dust, vibration, corrosive substances, explosive materials, high temperature and moisture.
- If the oximeter gets wet, please stop operating it.
- When it is carried from cold environment to warm or humid environment, please do not use it immediately.
- DO NOT operate keys on front panel with sharp materials.
- High temperature or high pressure steam disinfection of the oximeter is not permitted. Refer to User Manual in the relative chapter (7.1) for instructions of cleaning and disinfection.
- Do not submerge the oximeter in liquid. When it needs cleaning, please wipe its surface with medical alcohol by soft material. Do not spray any liquid on the device directly.
- When cleaning the device with water, the temperature should be lower than 140° F (60° C).

- Fingers which are too thin or too cold may affect the measurement accuracy. Please clip the thicker finger, such as thumb and middle finger deeply into the probe.
- The pulse oximeter can be used for adult and pediatric patients, based on the manufacturer supplied finger probe.
- The update period of data is less than 5 seconds, and is subject to change based on variations in individual pulse rates.
- Please read the measurements on the waveform screen once the readings have stabilized. This measure value is most accurate.
- If some abnormal conditions appear on the screen during test process, pull out the patient finger and reinsert to restore normal use. If continued abnormal conditions appear, contact your equipment provided.
- The device has a normal life expectancy of three years since the first electrified use.
- This device has alarm functions. Please reference chapter 6.1 for details.
- The device may not work for all patients. If you are unable to achieve stable readings, discontinue use.

2 OVERVIEW

The pulse oxygen saturation is the percentage of HbO_2 in the total Hb in the blood, so-called the O_2 concentration in the blood. It is an important bio-parameter for the respiration. A number of diseases relating to respiratory system may cause the decrease of SpO_2 in the blood, furthermore, some other causes such as the malfunction of human body's self-adjustment, damages during surgery, and the injuries caused by some medical checkup would also lead to the difficulty of oxygen supply in human body, and the corresponding symptoms would appear as a consequence, such as vertigo, impotence, vomit etc. Serious symptoms might bring danger to human's life. Therefore, prompt information of patients' SpO_2 is of great help for the doctor to discover the potential danger, and is of great importance in the clinical medical field.

The Handy-Ox Pulse Oximeter is a compact, efficient and convenient portable pulse oximeter. It is only necessary for a patient to put one finger into probe for diagnosis, and the display screen will directly show the measured value of SpO_2 and pulse rate with a high level of accuracy.

2.1 Features

- A. Operation of the product is simple and convenient.
- B. The product is compact, light weight and portable.
- C. Low power consumption.

The problem of overrating would emerge when the patient is suffering from toxicosis which caused by carbon monoxide, the device is not recommended to be used under this circumstance.

2.2 Environment Requirements

Storage Environment

- a) Temperature $-40^{\circ} \sim 140^{\circ} \text{ F}$ ($-40^{\circ} \sim 60^{\circ} \text{ C}$)
- b) Relative humidity: $\leq 95\%$
- c) Atmospheric pressure: $500\text{hPa} \sim 1060\text{hPa}$

Operating Environment

- a) Temperature: $50^{\circ} \sim 104^{\circ} \text{ F}$ ($10^{\circ} \sim 40^{\circ} \text{ C}$)
- b) Relative Humidity: $\leq 75\%$
- c) Atmospheric pressure: $700\text{hPa} \sim 1060\text{hPa}$

3. PRINCIPLE

Principle of the Oximeter is as follows: An experience formula of data process is established taking use of Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive Hemoglobin (Hb) and Oxyhemoglobin (HbO_2) in glow & near-infrared zones. Operation principle of the device is: Photoelectric Oxyhemoglobin Inspection Technology is adopted in accordance with Capacity Pulse Scanning & Recording Technology, so that two beams of different wavelength of lights can be focused onto human nail tip through perspective clamp finger-type sensor. Then measured signal can be obtained by a photosensitive element, information acquired through which will be shown on screen through treatment in electronic circuits and microprocessor.


Figure 1

4. TECHNICAL SPECIFICATIONS

4.1 Main Performance

- A. SpO₂ value display
- B. Pulse rate value display, bar graph display
- C. Pulse waveform display
- D. Low-voltage indication: low-voltage indicator will appear before operation and accuracy are affected.
- E. The display mode can be changed
- F. Screen brightness can be changed
- G. Pulse sound indication

H. Alarm function

I. With SpO₂ value and pulse rate value record function, the stored data can be uploaded to computer

J. It can be connected with an external oximeter probe

K. Real-time data can be transmitted to computer

4.2 Main Parameters

A. Measure of SpO₂

Measure range: 0%~100%

Accuracy:

When the SpO₂ measure range is 70%~100%,the permission of absolute error is $\pm 2\%$;
below 70% unspecified

B. Measure of pulse rate

Measure range:30bpm~250bpm

Accuracy: ± 2 bpm or $\pm 2\%$ (select larger)

C. Resolution

SpO₂: 1%, Pulse rate: 1bpm.

D. Measure Performance in Weak Filling Condition

SpO₂ and pulse rate can be shown correctly when pulse-filling ratio is 0.4%. SpO₂ error is $\pm 4\%$, pulse rate error is ± 2 bpm or $\pm 2\%$ (select larger).

E. Resistance to surrounding light

The deviation between the value measured in the condition of man-made light or indoor natural light and that of darkroom is less than $\pm 1\%$.

F. Power supply requirement: 3.6 V DC ~ 4.2V DC.

G. Optical Sensor

Red light (wavelength is 660nm, 6.65mW)

Infrared (wavelength is 880nm, 6.75mW)

H. Adjustable alarm range

SpO₂: 0%~100%

Pulse Rate: 0bpm~254bpm

5 INSTALLATION

5.1 View of the Front Panel


Figure 2
Front View

5.2 Underside View and Left View

Figure 3
Underside
and Left View

1. Probe jack: It is used to connect a fingertip probe to measure the oxygen saturation and pulse rate.
2. USB port: It is used to connect a personal computer to export the trend data(or real-time data) or charge the lithium battery.

5.4 Rear View

Figure 4
Rear View

(Refer to Figure 4 and insert the lithium battery properly in the right direction.)

5.4 Accessories

- A. A lithium battery
- B. A User Manual
- C. A power adapter
- D. USB Cord
- E. A disk (PC software)
- F. An adult-oximeter probe 18720-ADPROBE
An infant-oximeter probe (optional) 18720-PDPROBE

6 OPERATING GUIDE

6.1 Application Method

A. Measuring the data

- a) Open the lid of battery box and put the battery inside, then close.
- b) Put the fingertip probe into the jack on the underside of the oximeter. Only use Drive Medical supplied accessories.
- c) Put the finger into the probe.
- d) Press and hold the “power on/off button” for a few seconds. The device will power on within a few seconds and will start to self-check. The device displays the user interface after self-check.
- e) Do not shake the finger and keep the patient in a calm state during measurement.
- f) The data can be read directly from the screen in the user interface.

Fingernails and the luminescent tube should be at the same side.

If the alarm function is on, the device will provide a medium priority alarm signal when the probe plug or finger is out. An intermittent alarm will occur and the user interface presents “FINGER OUT”.

Medium priority indicates that prompt operator response is required.

The corner position of the plug should be inserted into the right place of the recessed portion of the socket, just as shown in the picture.

Figure 5

Figure 6

B. Change display direction

On the user interface, press the "screen change button" to change the display direction.

C. Pause alarm

- The Digital Handy-Ox will alarm if the PR or SpO₂ exceed preset low and high parameters, if the battery has a low charge, and if the patient's finger is removed.
- When alarm is on, press the "alarm pause button" to pause the alarm. In order to reset the alarm, press the "alarm pause button" again within 60s.
- If you want to turn off the alarm for good, enter the menu.

D. Menu operation

In the user interface, press the "menu button" to enter the menu options screen ((Figure 7). Users can adjust the setting through the main menu, such as alarm, pulse sound indication, PR and SpO₂ low and high alarm settings, backlight, data storage, and ID function.


Figure 7
Main Menu Interface

a) Alarm setting

In the main menu interface, press the "up button" or "down button", and move the menu choice bar to "Alarm", then press the "menu button" to enter the alarm setting menu of (Figure 8).

Alarm		
Spo2	ALH HI	99
Spo2	ALH LO	85
PR	ALH HI	150
PR	ALH LO	50
Alarm		ON
Exit		

Figure 8
Alarm Setting Menu

a. The high/low alarm limit setting

Press the "up button" or "down button" to choose the parameter to be adjusted, then press the "menu button" again to enter the similar dialog box shown as figure 9, then press the "up button" or "down button" to change data. Each press of the "up button" or "down button", the data will raise or descend in single increments. Continue until the desired setting is reached, then press the "menu button" to save the setting.


If the alarm function is on, the device will provide medium priority alarm signal when the data of SpO₂ or pulse rate is beyond the limit. Intermittent alarm will occur and the measure value becomes yellow.

Medium priority indicating that prompt operator response is required.

When setting up operation, the operator inserts the patient's finger, and sets the SpO2 ALM HI lower than the current measure value, then the machine will start alarm. The pulse rate alarm test is the same with the above.

Figure 9

b. The alarm state setting

Press the “up button” or “down button”, and move the menu choice bar to “Alarm”, then choose the alarm state “on/off” by pressing the “menu button”. Choose “on” to start the alarm function, and choose “off” to disable it.

b) Pulse sound indication setting

In the main menu interface, press the “up button” or “down button”, and move the menu choice bar to “Pulse Sound”, then choose the pulse sound indication state “on/off” by pressing the “menu button”. Choose “on” to start the function, and choose “off” to disable it.

c) Backlight adjustment

In the main menu interface, press the “up button” or “down button”, and move the menu choice bar to “Brightness”, then press the “menu button” to change the number in order to adjust the brightness of screen.

d) Data storage setting

This device has the function of 24 hours data storage, it has the capability of accurately storing pulse rate and SpO₂ data, then transmitting the data to a computer with the data line for review and analysis.

- a. In the main menu interface, press the “up button” or “down button”, and move the menu choice bar to “Record”, then press the “menu button” again to enter the dialog box of Figure 10. If it is in recording state, the “Stop” is displayed in the Menu; if it is not in recording state, the “Start” is displayed in the Menu.

Figure 10

- b. When the “Start” is displayed in the menu, if the memory still has old data in storage, the “Do you really want to recover the memory” dialog box will be displayed when pressing the “menu” button. Press the “up button” or “down button” to choose setting, and press the “menu button” to confirm setting, then the dialog box will be displayed as in Figure 11. To set the starting time for data storage by users, the user can change the setting value by pressing the “up button” or “down button”, and can confirm the setting value by pressing the “menu button”. The black cursor will move to next place after setting. When the black cursor moves to “Yes”(store) or “No”(don’t store), press the “up button” or “down button” to choose setting, then press the “menu button” to confirm setting.

Figure 11

- c. When the “Stop” is displayed in menu, the “Do you really want to stop recording” dialog box will be displayed when pressing the “menu button”. Press the “up button” or “down button” to choose setting, then press the “menu button” to confirm setting.
- d. If the data storage function is turned on, when return to the measure interface, a red “REC” sign and a flashing red dot would appear on screen, which means the device is in the state of storing.
- e. In the state of storing, whatever interface the device is in (measure interface, menu interface), the sign “Recording” would appear on the screen in 30 seconds, then the screen will be automatically shut down. If pressing any button (power on/off excluded) at this moment, the sign “Recording” would appear on the screen, then the screen will be automatically shut down again; if pressing the “power on/off button”, the device would return to the former interface.


- f. If turning on the data storage function, the old data will be automatically removed.
- g. In the state of data storing, after the screen is automatically shut down, the pulse sound indication would be off for saving power.
- h. When the storage space is full, it displays “Memory is full” on the screen, and then shut down in a few seconds. But it will still display “Memory is full” the next time you power on the device. If pressing any button (power on/off excluded) again, it will enter the measure interface.

e) **Stored data transmission setting**


- a. Connect the device to computer with data line, then open the PC software. Please refer to <SpO₂ Assistant user manual> for detailed operation method of uploading data. The software symbol is as Figure 12.

⚠ The computer the device is connected to should accord to the standard of IEC60950, and the system which is composed with the device accord with the requirements of IEC60601-1-1.

Figure 12
Software Symbol

- b. In the state of storing, it is not applicable for the users to upload the stored data to computer.
- c. When the stored data is being uploaded, the screen displays as figure 13.

Waiting...

Figure 13

d. When the uploading of stored data is finished, the device will automatically return to the measuring Interface.

f) Device ID

The user can set device ID by PC software. Please refer to <SpO₂ Assistant user manual> for details.

g) Exit the main menu

In the main menu interface, press the “up button” or “down button” in order to move the menu choice bar to “Exit”, then press the “menu button” to exit the main menu.

E. PC software operation

Please connect the device with computer by the USB data line, then double click “SpO₂ Assistant” icon to run the PC software. The functions such as uploading data and changing device ID could be carried out by the software. Please refer to <SpO₂ Assistant user manual> for detail.

If the user choose to turn on the synchronizing display function in computer, it would probably take several seconds for the data to appear on the computer screen. (If there is no data on the computer screen, unplug USB data line, then repeat step “E” again.)

F. Charge

There are two kinds of charge methods:

- Connect the device with computer by data line, then the device should be under charge state.
- Connect the device with power supply by power adaptor, then the device should be under charge state.
- The five status of battery power are shown as follows.

	Power supply by battery only, and battery status is full
	Battery status is not full
	
	Battery status is draining away
	Low power alarm indication (Please charge the battery)
the four status above shows dynamically in turn	Battery is charging

If the alarm function is on, the device will provide high priority alarm signal when the battery is in low power status. Intermittent alarm will occur and the battery icon turns red in the state of flashing. High priority indicating that immediate operator response is required.

6.2 Attention for Operation

- A. Please check the device before use, and confirm that it can work normally.
- B. The finger should be in a proper position (see the attached illustration of figure 6 for reference), or else it may result in inaccurate measure.
- C. The ray between luminescent tube and photoelectric receiving tube must get across subject's arteriole.
- D. The SpO₂ sensor should not be used at a location or limb tied with arterial canal or blood pressure cuff or receiving intravenous injection.
- E. Do not fix the SpO₂ sensor with adhesive, or else it may result in venous pulsation and inaccurate measure of SpO₂ and pulse rate.
- F. Excessive ambient light may affect the measure result. It includes fluorescent lamp, dual ruby light, infrared heater, direct sunlight and etc.
- G. Intense activity of the subject or extreme electrosurgical interference may also affect the accuracy.
- H. Patient should not use enamel or other makeup.
- I. Please clean and disinfect the device after operating according to the User Manual (7.1).

6.3 Clinical Restrictions

- A. As the measure is taken on the basis of arteriole pulse, substantial pulsating blood flow of subject is required. For a subject with weak pulse due to shock, low ambient/body temperature, major bleeding, or use of vascular contracting drug, the SpO₂ waveform (PLETH) will decrease. In this case, the measure will be more sensitive to interference.
- B. For those with a substantial amount of staining dilution drug (such as methylene blue, indigo green and acid indigo blue), or carbon monoxide hemoglobin (COHb), or methionine (Me+Hb) or thiosalicylic hemoglobin, and some with icterus problem, the SpO₂ determination by this device may be inaccurate.
- C. The drugs like dopamine, procaine, prilocaine, lidocaine and butacaine may also be a major factor resulted in serious error of SpO₂ measure.
- D. As the SpO₂ value serves as a reference value for judgement of anemic anoxia and toxic anoxia, and some patients with serious anemia may also report good SpO₂ measure value.

7 MAINTAIN, TRANSPORTATION AND STORAGE

7.1 Cleaning and Disinfecting

Use medical alcohol to disinfect the device. Dry or clean it with clean soft cloth.

7.2 Maintain

- A. Please clean and disinfect the device before use according to the User Manual (7.1).
- B. Please recharge the battery when the screen shows.
- C. Recharge the battery soon after the over-discharge. The device should be recharged every six months when it is not regularly used. It can extend the battery life following this guidance.
- D. Please take out the battery if the oximeter is not in use for a long time.
- E. Users are advised to calibrate the device annually, or in according to the calibration program of the hospital or the healthcare provider facility. It also can be performed at the state-appointed agent.

7.3 Transportation and Storage

- A. The packed device can be transported by ordinary conveyance. The device cannot be transported mixed with toxic, harmful, corrosive material.
- B. The packed device should be stored in room with no corrosive gases and good ventilation. Temperature: -40°F ~ 140°F(-40°C~60°C); Humidity: ≤95%.

8 TROUBLESHOOTING

Trouble	Possible Reason	Solution
The SpO ₂ and Pulse Rate cannot be displayed normally	<ol style="list-style-type: none">1. The finger is not properly positioned.2. The patient's SpO₂ is too low to be detected.	<ol style="list-style-type: none">1. Place the finger properly and try again.2. Try again; Go to a hospital for a diagnosis if you are sure the device is operating properly.
The SpO ₂ and Pulse Rate are not displayed stably	<ol style="list-style-type: none">1. The finger is not placed inside deep enough.2. The finger is shaking or the patient is moving.	<ol style="list-style-type: none">1. Place the finger properly and try again.2. Let the patient keep calm.
The device cannot be turned on	<ol style="list-style-type: none">1. The battery is drained away or almost drained away.2. The battery is installed incorrectly.3. The malfunction of the device.	<ol style="list-style-type: none">1. Please recharge the battery.2. Please Install the battery again.3. Please contact your equipment provider.
The display is off suddenly	The battery is drained away or almost drained away.	Please recharge the battery.
The device cannot be used for full time after charge	<ol style="list-style-type: none">1.The battery is not full charged.2. The battery is dead.	<ol style="list-style-type: none">1. Please recharge the battery.2. Please contact the local service center.
The battery cannot be full charged even after 10 hours charging time	The battery is dead.	Please contact the local service center.

9 KEY OF SYMBOLS

Signal	Description
	Refer to instruction manual/booklet
%SpO ₂	The pulse oxygen saturation(%)
PRbpm	Pulse rate (bpm)
	Close the alarm sound indication forever
	Pause the alarm sound indication
	Open the alarm sound indication
	Close the pulse sound indication
	Open the pulse sound indication
	Power on/off button
	Alarm pause button/up button

Signal	Description
	Menu button
	Screen change button/down button
	BF Type application part
SN	Serial number
	<ol style="list-style-type: none"> 1. The finger clip falls off (no finger inserted) 2. Probe error 3. Signal inadequacy indicator
IP22	International Protection
	WEEE (2002/96/EC)

10 FUCTION SPECIFICATION

Information	Display Mode
The Pulse Oxygen Saturation (SpO ₂)	2-digit digital TFT display
Pulse Rate (PR)	3-digit digital TFT display
Pulse Intensity (bar-graph)	Bar-graph TFT display
SpO₂ Parameter Specification	
Measuring range	0%~100%, (the resolution is 1%).
Accuracy	70%~100%:±2%,Below 70% unspecified.
Average value	Calculate the Average value in every 4 measure value. The deviation between average value and true value does not exceed 1%.
Pulse Parameter Specification	
Measuring range	30bpm~250bpm, (the resolution is 1bpm)
Accuracy	±2bpm or ±2% (select larger)
Average pulse rate	Moving calculate the Average pulse rate every 4 cardio-beat,s cycle. The deviation between average value and true value does not exceed 1%.
Safety Type	Interior Battery, BF Type application part
Pulse Intensity	
Range	Continuous bar-graph display, the higher display indicate the stronger pulse.

FUNCTION SPECIFICATION (CONTINUED)

Battery Requirement	
Voltage 3.7 rechargeable lithium battery × 1	
Battery working life	
Charge and discharge no less than 500 times.	
Power Adapter	
Input Voltage	100 to 240 VAC, 50/60 Hz
Output voltage	5 VDC
Output current	1A
Oximeter Probe	
Wavelength:660nm 880nm	
Dimensions and Weight	
Dimensions	8.7 (L) x 4.5 (W) x 2.2 (H) in (87(L) × 45(W) × 22(H) mm)
Weight	About 6.2 ounces (175 g) with the lithium battery

APPENDIX 1

Guidance and manufacturer's declaration electromagnetic emission for all EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration—electromagnetic emission		
The <i>Pulse Oximeter (18720)</i> is intended for use in the electromagnetic environment specified below. The customer or the user of the <i>Pulse Oximeter (18720)</i> should assure that it is used in such an environment		
Emission test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The <i>Pulse Oximeter (18720)</i> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The <i>Pulse Oximeter (18720)</i> is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emission IEC 61000-3-3	Not applicable	

Guidance and manufacturer's declaration—electromagnetic immunity for all EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration—electromagnetic emission			
The <i>Pulse Oximeter (18720)</i> is intended for use in the electromagnetic environment specified below. The customer or the user of the <i>Pulse Oximeter (18720)</i> should assure that it is used in such an environment			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6KV contact ±8KV air	±6KV contact ±8KV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment

Guidance and manufacturer's declaration-electromagnetic immunity for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacturer's declaration-electromagnetic immunity			
The <i>Pulse Oximeter (18720)</i> is intended for use in the electromagnetic environment specified below. The customer or the user of the <i>Pulse Oximeter (18720)</i> should assure that it is used in such an environment			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment -guidance
Conducted RF IEC 61000-4-6	3Vrms 150kHz to 80MHZ	3Vrms	<p>Portable and mobile RF communication equipment should be used no closer to any part of the <i>Pulse Oximeter (18720)</i>, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_i} \right] \sqrt{P} \quad 80\text{MHz to } 800\text{MHz}$ $d = \left[\frac{7}{E_i} \right] \sqrt{P} \quad 800\text{MHz to } 2.5\text{GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>
Radiated RF IEC 61000-4-3	3V/m 80MHz to 2.5GHz	3V/m	
NOTE 1 At 80MHz and 800MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
<p>a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which The <i>Pulse Oximeter (18720)</i> is used exceeds the applicable RF compliance level above, the <i>Pulse Oximeter (18720)</i> should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the <i>Pulse Oximeter (18720)</i>.</p> <p>b. Over the frequency range 150 KHz to 80 MHz, field strengths should be less than 3V/m.</p>			

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM for EQUIPMENT or SYSTEM that not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the Pulse Oximeter (18720)		
The Pulse Oximeter (18720) is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Pulse Oximeter (18720) can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Pulse Oximeter (18720) as recommended below, according to the maximum output power of the communications equipment.		
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)	
	80MHz to 800MHz $d = \left[\frac{3.5}{E_f} \right] \sqrt{P}$	800MHz to 2.5GHz $d = \left[\frac{7}{E_f} \right] \sqrt{P}$
0.01	0.1167	0.2334
0.1	0.3689	0.7378
1	1.1667	2.3334
10	3.6893	7.3786
100	11.6667	23.3334
For transmitters rated at a maximum output power not listed above, the recommended separation distance in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80MHz and 800MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		

APPENDIX 2

State	Alarm condition delay	Alarm signal generation delay
Low voltage alarm	1s	20ms
SpO ₂ alarm	330ms	20ms
Pulse rate alarm	330ms	20ms
Probe error alarm	16ms	20ms

WARRANTY

Your Drive brand product is warranted to be free of defects in materials and workmanship for two years from the original purchase date. The device was built to exacting standards and carefully inspected prior to shipment. This 2 year Limited Liability warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service. In the event of a defect covered by this warranty, we will, at our option, repair or replace the device.

This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. If you have questions about your Drive device, or this warranty, please contact an authorized Drive Medical provider.

Drive Medical Design & Manufacturing

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