**WelchAllyn**

9600 Plus Calibration Tester

Copyright 2004, Welch Allyn, Inc. All rights are reserved. No one is permitted to reproduce or duplicate, in any form, this manual or any part thereof without permission from Welch Allyn.

Welch Allyn assumes no responsibility for any injury to anyone, or for any illegal or improper use of the product, that may result from failure to use this product in accordance with the instructions, cautions, warnings, or statement of intended use published in this manual.

Welch Allyn is a registered trademark of Welch Allyn. Braun ThermoScan is a registered trademarks of Braun Corporation.

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## Introduction

---

The 9600 Plus Calibration Tester verifies calibration of the following thermometers:

- Welch Allyn SureTemp<sup>®</sup> models (678, 678 OEM, 679 and 986)
- Welch Allyn SureTemp<sup>®</sup> Plus models (690, 692 and 692 OEM)
- Braun ThermoScan<sup>®</sup> Pro 3000 and Pro 4000

The 9600 Plus Calibration Tester provides a tightly regulated temperature source that eliminates the setup, precision calibration, and maintenance requirements of water baths and other general purpose calibration equipment. The user takes the “temperature” of the 9600 Plus Calibration Tester with the knowledge and assurance that precise readings are obtained because the unit is factory calibrated to National Institute of Standards and Technology (NIST) or equivalent devices. User adjustments are not required or recommended.

The “blackbody” of the 9600 Plus Calibration System acts as a heat reservoir for the thermistor probe and as a heat radiator for the infrared ear thermometer.

The unit contains solid-state circuitry for long-term stability with minimal maintenance. Welch Allyn uses only precision components in the construction of the 9600 Plus Calibration Tester. To maintain traceability to NIST or equivalent, return the 9600 Plus Calibration Tester to the factory for annual calibration.

## Symbols and Descriptions

Double Insulation.

Direct Current.

Handle with Care.

Conforms to IEC 61010-1:2001 and EN 61010-1:2001

Attention! Caution; Consult Directions for Use for more information.

Transport/Storage Temperature

Relative Humidity

Separate collection for electronic equipment

Recyclable

## Cautions

A caution statement in this manual identifies a condition or practice which, if not corrected or discontinued immediately, could lead to equipment failure, equipment damage, or data loss.

**Caution** Do not immerse or soak the 9600 Plus Calibration Tester in any type of fluid.

**Caution** Only trained Biomedical personnel should operate the 9600 Plus Calibration Tester and conduct calibration verification of thermometers.

**Caution** When not using the Thermistor or Ear Device Ports, keep them covered with the Dust Cover.

**Caution** Store thermometers for testing in the same room as the 9600 Plus Calibration Tester for approximately 30 minutes prior to testing to allow for thermal accommodation.

**Caution** Do not clean the Thermistor or Ear Device Ports. Doing so may damage the surface. Only a Welch Allyn authorized Service Center should clean the Device Ports.

## Unpacking the 9600 Plus Calibration Tester

Unpack the 9600 Plus Calibration Tester and inspect for missing items. Retain the shipping materials to return the device to Welch Allyn for recalibration.

The following items are enclosed:


- A 9600 Plus Calibration Tester
- A package containing a power supply and four mains plug inserts
- The Directions for Use
- A Carrying Case

If any of these items are missing after opening, please contact Welch Allyn Customer Service at the number provided on [page 28](#).

## 2 Controls and Indicators

### LCD Display

Figure 1. LCD Display and Front Buttons

### Calibration Set Point Temperature

There are three Calibration Set Point temperatures used to test the calibration of the thermometer. The three temperatures are 36° C (96.8° F), 38.5° C (101.3° F), and 41°

C (105.8° F). The LCD displays the current Calibration Set Point Temperature in the upper left corner of the LCD.

## Current Cavity Temperature

The current temperature of the 9600 Plus Calibration Tester is indicated by the large numbers in the center of the display (Figure 1).

As the device is transitioning to the current Calibration Set Point Temperature, the actual temperature flashes on the display. When the device reaches the Calibration Set Point Temperature, the temperature in the center of the display stops flashing.

## Temperature Units

The 9600 Plus Calibration Tester is capable of displaying the temperature in degrees Fahrenheit or Celsius.

## Set Point Mode

The LCD will display either a "D" (Default) or "L" (Last Used) in the upper left hand corner (Figure 1). The 9600 Plus Calibration Tester has two set point modes of operation. The Default mode (from the factory) causes the 9600 Plus Calibration Tester to heat up to the lowest set point (36.0°C / 96.8°F) upon startup. The Last Used mode causes the 9600 Plus Calibration Tester to heat up to the last temperature setting that was used.

To change the Set Point Mode, press and hold the C/F and Arrow buttons below the LCD display simultaneously until a single audio beep is heard. A "D" or an "L" will appear in the upper left corner of the display for the new mode.

## Power On Sequence

The 9600 Plus Calibration Tester receives power from a wall-mount power pack that plugs into the back of the unit. There is no ON/OFF switch. On power up, the 9600 Plus Calibration Tester performs internal self-checks, displays an all-segments test pattern, beeps and displays the device's current temperature. The 9600 Plus Calibration Tester immediately applies heat, if needed, to reach the selected Calibration Set Point Temperature.

## Front of the 9600 Plus Calibration Tester

### Temperature Selection Button


This button allows the user to select one of three set point temperatures: 36.0 °C (96.8 °F), 38.5 °C (101.3 °F), or 41.0 °C (105.8 °F).

To scroll from one set point to the next, press and hold the Temperature Selection button until a beep is heard. The newly selected set point appears in the upper left corner of the LCD display. The device's current temperature is displayed, will start to flash, and continue flashing until the cavity reaches the equilibrium at the new set point.

### Celsius/Fahrenheit Selection Button

This button allows the user to select degrees Celsius or degrees Fahrenheit. To change the temperature units, press and hold the C/F button until a beep is heard. The temperatures will appear in either "°C" or "°F"; see [\(Figure 1\)](#).

**Figure 2. Front of the 9600 Plus Calibration Tester**

## Dust Cover

This black, plastic cover keeps dust from entering the Thermistor and Ear Device Ports. To expose the Thermistor and Ear Device Ports, rotate the Dust Cover counter clockwise. To cover the ports, rotate the Dust Cover clockwise.

**Caution** When not using the Thermistor or Ear Device Ports, keep them covered with the Dust Cover.

## Thermistor and Ear Device Ports

The Ear Device Port is the larger of the two ports located under the dust cover and is for the Braun ThermoScan Pro 3000 and Pro 4000 ear thermometers ([Figure 2](#)).

The Thermistor Device Port is the smaller port located under the dust cover and is for the SureTemp and SureTemp Plus thermometers.

The surface texture inside the Ear Device Port is painted with a special dull, black paint to enable Ear device readings. You can only see its surface under bright light. The Thermistor Device Port is not painted black. Keep both clean and dry.

The coating in the device ports of the unit is fragile, so please follow the cleaning instructions carefully. Do not touch the surface, and do not allow any foreign objects to be inserted into the ports. Avoid spilling liquids onto the unit. If liquid or other contaminants enter the unit, Welch Allyn can not assure proper calibration and you must return the unit to Welch Allyn for recalibration.

## **3** Operation

---

### 9600 Plus Calibration Tester Setup

Place the 9600 Plus Calibration Tester on a level surface away from sunlight, drafts, and other sources of heat or cold. The 9600 Plus Calibration Tester takes approximately 20 minutes to heat to the lowest setting. When testing several thermometers at all three temperatures, it is recommended to test all probes at one Calibration Set Point Temperature before proceeding to the next Calibration Set Point Temperature. To further expedite testing start at the lowest Calibration Set Point Temperature. The 9600 Plus Calibration Tester does not have an internal fan, causing the cool down time to be greater than the time it takes to heat up.

**Figure 3. Power Adapter and Mains Plug Inserts**

For the setup:


1. Choose the proper mains plug insert and slide it over the two prongs in the power converter (Figure 3).
2. Plug the power adapter into the 9600 Plus Calibration Tester (Figure 2) and the opposite end into a wall outlet.
3. Observe the Set Point Mode in the upper left hand corner of the LCD display (Figure 1). If the unit displays a "D", it is in Default Mode and will heat to the lowest Set Point Temperature. If you do not want to conduct testing at this Set Point Temperature, press and hold the Temperature Selection button to select the desired setting. The temperature display will flash before staying on continuously to indicate the 9600 Plus Calibration Tester has stabilized and is ready for use.

## Welch Allyn SureTemp and SureTemp Plus Thermometers

**Caution** Store thermometers for testing in the same room as the 9600 Plus Calibration Tester for approximately 30 minutes prior to testing to allow for thermal accommodation.

**Figure 4. 9600 Plus Calibration Tester with the Welch Allyn SureTemp Plus Thermometer**

To begin calibration verification of the thermometers:


1. Remove the probe from the probe well and clean it with either a 70% isopropyl alcohol solution, a 10% chlorine bleach solution, or a non-staining disinfectant. Let the probe air dry. Do not apply a probe cover.
2. Place the thermometer in Monitor Mode, refer to the thermometer's Operator's Manual.
3. Insert the probe into the Thermistor Device Port ([Figure 4](#)).
4. Wait for approximately one minute or until temperature on the thermometer is stable for ten seconds. Compare the thermometer's temperature reading to the 9600 Plus' Calibration Set Point Temperature. If the temperatures are within  $\pm 0.1^{\circ}\text{C}$  ( $\pm 0.2^{\circ}\text{F}$ ), the thermometer is within calibration.
5. Test all available thermometers for calibration verification at the current Calibration Set Point Temperature. Proceed to the next Calibration Set Point Temperature, see "[Changing the Calibration Set Point Temperature](#)" ([page 19](#)).

## Braun ThermoScan Pro 4000 Model

**Caution** Store thermometers for testing in the same room as the 9600 Plus Calibration Tester for approximately 30 minutes prior to testing to allow for thermal accommodation.

**Figure 5. 9600 Plus Calibration Tester with the Braun ThermoScan Pro 4000**

To begin calibration verification of the thermometers:

1. Clean the probe tip with a cotton swab slightly moistened with isopropyl alcohol, remove excess alcohol with a clean

cotton swab, and let air dry for 5 minutes. Do not use any chemical other than alcohol to clean the probe window.


2. Place the Braun ThermoScan Pro 4000 in Calibration Mode using the following steps.
  - a. Push and release the I/O MEM button to turn the product on. The display will show symbols and functions. The Pro 4000 will perform an automatic self check. After a few seconds you will hear a beep and see three dashed lines, a sound symbol, and ° C or ° F on the display.
  - b. Push the I/O MEM button again and keep it pressed for the following steps:
    - After approx. 3 seconds the "OFF" symbol will flash on the display (keep the button depressed).
    - When you hear a beep, release the button immediately.
    - The Pro 4000 is now in calibration check mode and the display is flashing and showing the "CAL" symbol.
3. Apply a new probe cover. Place the probe firmly into the Ear Device Port ([Figure 5](#)).
4. Wait approximately three seconds, press the thermometer Start button, and watch for the "Exac Temp" light to flash.
5. Leave the thermometer in the 9600 Plus Calibration Tester until you hear a beep.
6. Remove the Pro 4000 Thermometer from the 9600 Plus Calibration Tester and read the temperature in the thermometer's display. If the temperatures are within  $\pm 0.2^{\circ}\text{C}$  ( $\pm 0.4^{\circ}\text{F}$ ), the thermometer is within calibration.

7. Wait one full minute before taking another reading with the same thermometer. Repeated measurements in short sequence may cause higher readings.
8. The device will exit CAL mode after four minutes.
9. Test all available thermometers for calibration verification at the current Calibration Set Point Temperature. Proceed to the next Calibration Set Point Temperature, see ["Changing the Calibration Set Point Temperature"](#) (page 19).

## Braun ThermoScan Pro 3000 Model

**Caution** Store thermometers for testing in the same room as the 9600 Plus Calibration Tester for approximately 30 minutes prior to testing to allow for thermal accommodation.

**Figure 6. 9600 Plus Calibration Tester with Braun ThermoScan Pro 3000**

1. Clean the probe tip with a cotton swab slightly moistened with isopropyl alcohol, remove excess alcohol with a clean cotton swab, and let air dry for 5 minutes. Do not use any chemical other than alcohol to clean the probe window.
2. Place the Braun ThermoScan Pro 3000 in calibration check mode using the following steps.
  - a. Wait until the thermometer is in Standby mode (the screen is blank).
  - b. Press down and hold the Activation button, do not release until [Step d.](#)
  - c. Press the  /Mem button.

- d. The LCD will flash "CAL" after the self check sequence. The display will flash "CAL" even if you have taken a temperature.
  - e. Release both buttons.
3. Apply a new probe cover. Place the device firmly into the Ear Device Port ([Figure 6](#)).
  4. Wait approximately three seconds, press the thermometer Activation button, and listen for the beep.
  5. Remove the Pro 3000 Thermometer from the 9600 Plus Calibration Tester and read the temperature in the thermometer's display. If the temperatures are within  $\pm 0.2^{\circ}\text{C}$  ( $\pm 0.4^{\circ}\text{F}$ ), the thermometer is within calibration.
  6. Wait one full minute before taking another reading with the same thermometer. Repeated measurements in short sequence may cause higher readings.
  7. The device will exit Cal mode after approximately 4 minutes.
  8. Test all available thermometers for calibration verification at the current Calibration Set Point Temperature. Proceed to the next Calibration Set Point Temperature, see "[Changing the Calibration Set Point Temperature](#)" ([page 19](#)).

## Changing the Calibration Set Point Temperature

To scroll from one set point to the next, press and hold the Temperature Selection button until a beep is heard. The newly selected set point appears in the upper left corner of the LCD display. The device's current temperature is displayed, will start to flash, and continue flashing until the cavity reaches the equilibrium at the new set point.

## 4 Maintenance

---

### Cleaning

To clean the 9600 Plus Calibration Tester:

1. Unplug the electrical cord from the back of the unit.
2. Using a cloth slightly dampened with a solution of warm water and mild soap, clean the outside of the unit. Take extra care so that soap and water do not enter the two device ports.
3. Rinse the cloth with water and thoroughly wring it out.
4. Wipe off the soap and water from the 9600 Plus.
5. Let the 9600 Plus Calibration Tester dry for 15 minutes.
6. Plug the electrical cord into the back of the unit.

**Caution** Do not clean the Thermistor or Ear Device Ports. Doing so may damage the surface. Only a Welch Allyn authorized Service Center should clean the Device Ports.

## Problem Solving

<b>Error</b>	<b>Problem</b>	<b>Solution</b>
E ▲	Calibration Tester cannot reach Set Point Temperature because room or device is too warm.	Ambient temperature and/or device temperature should be within 18.3°C to 26.7°C (65°F to 80°F). Allow 30 minutes for thermal accommodation prior to testing.
E ▼	Calibration Tester cannot reach Set Point Temperature because room or device is too cold.	Ambient temperature and/or device temperature should be within 18.3°C to 26.7°C (65°F to 80°F). Allow 30 minutes for thermal accommodation prior to testing.
Voltage	Power adapter voltage problem	Call Welch Allyn for replacement
E0.1	Broken Heater, not heating at all	Return to Welch Allyn for service
E0.2	Heater run away (heater does not turn off at appropriate time)	Return to Welch Allyn for service
E0.3	Open thermistor	Return to Welch Allyn for service
E0.4	Shorted thermistor	Return to Welch Allyn for service
E0.5	Open RCAL resistor	Return to Welch Allyn for service
E0.6	Shorted RCAL resistor	Return to Welch Allyn for service
E0.7	Open PTB resistor	Return to Welch Allyn for service
E0.8	Shorted PTB resistor	Return to Welch Allyn for service
E0.9	Catastrophic internal failure	Return to Welch Allyn for service

<b>Error</b>	<b>Problem</b>	<b>Solution</b>
E1.1	Device not calibrated	Return to Welch Allyn for service
E1.3	EEPROM failure	Return to Welch Allyn for service

## 9600 Plus Calibration Tester Checks

If several thermometers do not pass the calibration verification tests, check the 9600 Plus Calibration Tester to determine if it is out of calibration. Before sending the unit to an authorized Welch Allyn service center, complete the checks listed below.

1. Check the calibration date on the side of the 9600 Plus Calibration Tester. If the 9600 Plus Calibration Tester was calibrated over 12 months ago, do not use it. Send it to an authorized Welch Allyn Service center immediately for calibration.
2. Ensure no one has tampered with the 9600 Plus Calibration Tester. If someone has tampered with the unit or the calibration label(s), the calibration settings and warranty are void. Send it to an authorized Welch Allyn Service center immediately for calibration.
3. Ensure the 9600 Plus Calibration Tester is not exposed to extreme temperatures that can damage the precision thermistor. Exposure to high heat can cause the resistance versus temperature characteristic of the unit to shift, causing the 9600 Plus Calibration Tester to become inaccurate. Send it to an authorized Welch Allyn Service center immediately for calibration.
4. Verify the probe tip of the Braun ThermoScan Pro 3000 and the Braun ThermoScan Pro 4000 are clean. Always use new probe covers when calibration testing these thermometers.

5. Verify the ear thermometers are in calibration mode.
6. Verify the thermistor thermometers are in monitor mode.

## Environmental Checks

1. Verify that there is no excessive moisture or steam in the air. High humidity (above 90%) can cause electrical leakages that can temporarily shift the device out of calibration. If the thermometers and/or the 9600 Plus Calibration Tester are exposed to such conditions, allow each device to sit at normal room temperature and humidity for at least 30 minutes before conducting a calibration test.
2. Avoid exposure to extreme temperatures. Thermal gradients in the thermometer or the 9600 Plus Calibration Tester can cause unpredictable effects. Verify the 9600 Plus Calibration Tester is within the specified ambient temperature range of 18.3° C (65.0° F) to 26.7° C (80.0° F).

## Specifications

<b>Power Requirements</b>	Voltage: 100 - 240 V 50/60 Hz
<b>Current</b>	Heating: 480 mA After Warm-Up: 185 mA
<b>Power Consumption</b>	Heating: 2.4 W After Warm-Up: 0.92 W
<b>Calibration Tester</b>	Length: 254 mm (10.00") Width: 105 mm (4.13") Height : 85 mm (3.34") Weight (without Power Pack): 0.75 kg (1.65 Lbs)
<b>Power Pack (approximate dimensions)</b>	Length: 86 mm (3.38") Width: 51 mm (2.00") Height (with plug inserted): 62 mm (2.44") Weight: 0.18 kg (0.40 lbs)
<b>Environment</b>	Storage Temperature: 0°C – 50°C (32°F – 122°F) Ambient Operating Temperature: 18.3° C to 26.7°C (65° F to 80° F) Operating Humidity: 20% to 95% RH. Non-Condensing
<b>Operating Specifications</b>	Controlled Temperature of 9600 Plus Low Setting: 36°C (96.8°F) Intermediate Setting: 38.5°C (101.3°F) High Setting: 41°C (105.8°F) Temperature Accuracy: ±0.05° C (0.1°F) Maximum Pressure/Altitude: ≤ 2000 meters.
<b>Heating and Cooling Times</b> Tested at an ambient of 23.9°C (75°F)	Room Temperature to Low Setting: Approximately 20 min. Low Setting to Intermediate Setting : Approximately 6 min. Intermediate Setting to High Setting: Approximately 6 min. Cool from High Setting to Low Setting: Approximately 25 min.
<b>Installation Category II Pollution Degree II</b>	

## Electro-Magnetic Compatibility

This product has been tested and found to comply with all applicable requirements of council directive 89/336/EEC and its amendments.

This product has been tested and found to comply with FCC Part 15.107, Class B. These limits are designed to provide reasonable protection against harmful interference in a typical (medical) installation. This equipment generates, uses and can radiate radio frequency energy and if not used in accordance with the instructions, may cause harmful interference to other devices in close proximity. However, there is no guarantee that interference will not occur in a particular setting. If this equipment does cause harmful interference with other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device
- Increase the separation between equipment
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.
- Consult the manufacturer or field service technician for help.

## **5** Warranty and Service

---

### Limited Warranty

Instrumentation purchased new from Welch Allyn is warranted to be free from defects in material and workmanship under normal use and service for a period of one year from the date of first shipment from Welch Allyn. This warranty shall be fulfilled by Welch Allyn or its authorized representative repairing or replacing, at Welch Allyn's discretion, any such defect free of charge for parts and labor.

Welch Allyn should be notified via telephone of any defective product and the item should be immediately returned, securely packaged and postage prepaid to Welch Allyn. Loss or damage in shipment shall be at purchasers' risk.

Welch Allyn will not be responsible for loss associated with the use of any Welch Allyn product that (1) has had the serial number defaced, (2) has been repaired by anyone other than an authorized Welch Allyn Service Representative, (3) has been altered, or (4) has been used in any manner other than in accordance with instructions.

THIS WARRANTY IS EXCLUSIVE AND IN LIEU OF ANY IMPLIED WARRANTY OF MERCHANTABILITY, FITNESS FOR PARTICULAR PURPOSE, OR OTHER WARRANTY OF QUALITY, WHETHER EXPRESS OR IMPLIED, Welch Allyn WILL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES.

## Service

Before returning a product for repair or annual calibration you must obtain authorization from Welch Allyn. An RMA (Return Merchandise Authorization) number is given to you from our service personnel. Be sure to note this number on the outside of your shipping box. Returns without an RMA number will not be accepted for delivery. Please contact your customer service representative at the appropriate Welch Allyn service center listed below for further information.

All equipment used in calibrating the 9600 Plus Calibration Tester is traceable to NIST.

Welch Allyn assumes no responsibility if the 9600 Plus Calibration Tester is calibrated at an unauthorized service center. Send the 9600 Plus Calibration Tester to any one of the following authorized service centers for recalibration.

**USA, Central America, and South America**

Welch Allyn, Inc.  
4341 State Street Road  
Skaneateles Falls, NY 13153-0220  
Tel.: 1-800-535-6663  
Fax: 1-315-685-3361

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**Welch Allyn GmbH**