

## DeVilbiss Traveler® Portable Compressor Nebulizer System

### The **perfect** travel companion



DeVilbiss Traveler Portable  
Compressor Nebulizer Unit  
with battery attached  
Model 6910P-DR

The **DeVilbiss Traveler Portable Compressor Nebulizer System** is designed to help active customers get more out of their aerosol therapy, and **more out of life**. It's never been easier or more convenient to nebulize outside of the home.

The DeVilbiss Traveler includes compressor unit, **disposable** and reusable nebulizers, **five extra filters**, tubing adapter, adult and **pediatric** masks, AC adapter, DC adapter, carrying case, instruction guide and **instructional DVD**, as well as a **full-year** battery warranty. Units can be ordered with **optional battery**.

### Performance Features

- Delivers a **consistent and concentrated treatment** of small particles that penetrate deeply into the lungs, delivering medication where it's needed most
- **Fast treatment time** (approximately 10 minutes) with DeVilbiss 900 Series Reusable Nebulizer encourages patient compliance
- Can be used with **any manufacturer's nebulizer** – treatment times will vary
- **45-minute battery life** far outlasts treatment time

