



Height Measuring Rod

The AdirMed 970-01 Telescopic Height measuring Rod is a classic favorite. The headpiece of the 970-01 Telescopic Stadiometer is secured on both sides for accuracy and folds down for a greater degree of safety. In addition, the 970-01 Height Rod is double-graduated in inches and centimeters. When the 970-01 Telescopic Measuring Rod is added to a body scale, the height rod becomes a precision weighing and measuring package.

Standard Features

- To be used in 3 position:
- For measuring standing on weighing platform
- For measuring standing on the floor
- For measuring standing at the wall
- Measuring unit: cm + inch

Specifications

Measuring range: 600-2100 mm

Graduation: 1 mm

Dimensions: 93 x 35 x 35 mm

Weight approx.: 750 g

