



VIVA Slim Line Exam Table

ADI996-FBASE-M2-BLK

ADI996-FBASE-M2-BLU



Viva Slim Line Exam Table, Black

Our Viva Slim Line Exam Tables' compact design gives you that extra space you need for everyone to feel comfortable. The base of our tables has a dual-sided, connecting cabinet, allowing you to position the table however you like without compromising on cabinet entry. Two euro-style drawers with stainless steel slides located in the front of the base provide ample storage for your frequently used items. We've designed the drawers and door of the tables to be soft closing for that added touch of comfort. A pull-out step stool is also located in front of the table with a safety tread so your patients can safely get on and off the exam table.

We wanted to deliver a product you could feel confident will last; that is why the Viva Slim Line Exam Table frame is laminated and has chip-resistant edges, fortifying its longevity. Our tables have a built-in paper dispenser compatible with 17.7" – 20.8" paper rolls, so you can quickly change the liner between patients. Our provided antimicrobial mattress is designed to provide a safe and hygienic environment for your patients by reducing the risk of cross-contamination. Our Slim Line Exam Tables extensive features make them the ideal medical tables you can trust to give you ample storage while not compromising on space offering everything your patients need.

FEATURES

Accommodating: Adjustable backrest that can be positioned at different heights.

Build: The all-laminated frame maintains the longevity of the exam bed.

Convenience: A pull-out step stool is provided so patients can effortlessly get on and off the exam table.

Built-In: The included paper roll dispenser and paper cutter allows for quick changing of the liner.

Storage: Dual-sided cabinet and two front-facing euro-style drawers with stainless steel slides ensure plenty of space.

Sanitation: The antimicrobial mattress is designed to provide a safe and hygienic environment for patients by reducing the risk of cross-contamination.

SPECIFICATIONS

Base Weight	197 lb.
Base Dimensions	29" H x 47" W x 24.4" D
Cabinet Dimensions	19.7" H x 21.5" W
Material	Metal and Melamine Finished Board
Finish	Laminate Board
Length with Leg Rest Extended	72.04"
Leg Rest	15.7"
Weight Capacity	500 lb.
Compatible Paper Roll Size	Length: 17.71" – 20.86", Height: 2.99"
Step Stool Dimensions	16" D x 21.3" W x 9.8" L
Ground To Top of Step Stool	9.8"
Certification	California CARB II Compliant
Mattress Weight	35.2 lb.
Mattress	2.75-inch Padded Vinyl Antimicrobial Upholstery
Mattress Dimensions	2.75" x 58" x 27"
Total Weight	232.2 lb.
Color Options	Black / Midnight Express

FEATURES



Two Drawers




Retractable Step Stool



Dual-Sided Cabinet



Pull-out Leg Rest



Paper Roll Storage



Paper Roll Dispenser and Cutter

ADI996-FBASE-M2-BLU




**Viva Slim Line Exam Table,
Midnight Express**




Adjustable Back Rest

UPHOLSTERY COLOR OPTIONS



Black



*Midnight
Express*