

Infrared Thermometry



Part Number: TT-D-IR

Features:

- Temperatures in under 2 seconds
- Hygienic touchless loading and ejecting of probe covers
- Battery life of 15,000 measurements
- Intuitive icon based interface
- Pulse timer
- Tapered tip design suitable for babies, children and adults
- Durable stainless steel probe tip
- Convenient probe storage in base

Our infrared thermometer provides a quick and accurate non-invasive measurement. The tapered probe design and under 2 second measurement make it easy to use with babies and children.

3YEAR
warranty



Infrared Thermometry LCD Screen

Infrared Thermometry Options

Disposable Probe Covers



- Single use probe covers minimize the risk of cross contamination.
- Probe covers conveniently fit into holder.
- Touchless probe cover loading and ejecting support infection control.

Part Number: 303030: Pack of 2,112

Mounting Options



- Thermometry unit can be mounted on wall, rail system or roll stand.
- Mounting solution can be fully customized to your needs.

Specifications

Accuracy

- Quick Mode ± 0.1 °C (± 0.6 °F)
- Standard Mode ± 0.1 °C (± 0.2 °F)
- Direct Mode ± 0.1 °C (± 0.2 °F)

Measurement Response Time

- Quick Mode: 1-2 sec

Ambient Operating Range

- 16°C-33°C (60.8°F-107.6°F) at 10%-95% RH

Temperature Range

- 33.0°C-42.0°C (91.4°F-107.6°F)

Storage Temperature Range

- -25°C-55°C (-13°F-131°F) at 95% RH

Battery Life

- 15,000 measurements on 3 AAA batteries

Anti-Theft System

- Wall mount bracket

Infrared Thermometry Part Matrix

TT-X-IR

Mounting Options

D

Wall/Desktop Mount

R

Roll Stand


A

Amico Horizontal Rail Mount

V

Vertical Rail Mount

Distributed By:



AD-INF-THRM PRINTED JUL 13

Amico Diagnostic Incorporated 85 Fulton Way, Richmond Hill, ON L4B 2N4, Canada

Toll Free Tel: 1.877.264.2697 | Tel: 905.763.7778 | Fax: 905.763.8587

www.amico.com