

ELEVATES 75°



*Also available in 40"W Model:  
AM-BA240 (No Face Opening) and  
27"W Model: AM-BA227 (No Face  
Opening) with 400 lb. Capacity*

## AM-BA234 Bariatric Treatment Table

- 2-Section Bariatric Treatment Table
- 34"W x 76"L
- Min. Height: 18"
- Max Height: 37"
- 1 3/4" Firm Density Foam
- Backrest: 27" x 34" with Contoured Face Opening
- Body: 49" x 34"
- 500 lb. Capacity
- 4 Total Lock Swivel Casters

ELEVATES 78°



RAISES 45°  
LOWERS 90°

## AM-BA334 Bariatric Treatment Table

*Also available in 40"W  
Model: AM-BA340*

- 3-Section Bariatric Treatment Table
- 34"W x 76"L
- Min. Height: 18"
- Max Height: 37"
- 1 3/4" Firm Density Foam
- Head: 15" x 34" with Contoured Face Opening
- Center: 22" x 34"
- Foot: 39" x 34"
- 500 lb. Capacity
- 4 Total Lock Swivel Casters