

**P R O T E C T E D .**

**SUPERIOR FLUSHING,  
IMPROVED ACCESS**

**THE B-D® BLUNT PLASTIC CANNULA**



**Interlink®**

**MAKE THE CONNECTION**

**BECTON  
DICKINSON**


# SETTING THE STANDARD FOR IV ACCESS PERFORMANCE

Available separately and pre-attached to BD syringes

Clear material for enhanced visualization

Dual side ports for greater fluid dispersion in flushing and drug delivery

Center-point design with tapered tip for less resistance during insertion



THE NEW B-D® BLUNT PLASTIC CANNULA FEATURES A DYNAMIC DESIGN TO IMPROVE THE PERFORMANCE OF THE INTERLINK® IV ACCESS SYSTEM WHILE ELIMINATING THE NEED FOR SHARP NEEDLES.

THE CANNULA FEATURES DUAL SIDE PORTS WHICH CREATE GREATER FLUID DISPERSION IN FLUSHING AND DRUG DELIVERY. THE GREATER FLUID DISPERSION HELPS MINIMIZE THE POTENTIAL RISK OF INFECTION BY REDUCING THE RELIANCE ON TURBULENT FLUSHING TECHNIQUES.

THE CANNULA ALSO HAS A PROPRIETARY CENTER-POINT DESIGN TO PROVIDE SIMPLIFIED ACCESS TO THE INTERLINK® INJECTION SITE WITH SMOOTHER INSERTION AND REDUCED PENETRATION FORCE.

IT HAS CLINICAL EQUIVALENCE TO OTHER B-D CANNULA PRODUCTS FOR IV ACCESS APPLICATIONS, INCLUDING BLOOD COLLECTION.

THE INTERLINK® SYSTEM WAS THE FIRST TECHNOLOGY TO ALLOW FOR NEEDLELESS IV ACCESS. THE B-D® BLUNT PLASTIC CANNULA IS THE LATEST DEVELOPMENT IN OUR FAMILY OF PRODUCTS DESIGNED FOR BETTER INFUSION THERAPY.



B-D® BLUNT PLASTIC CANNULA ORDER INFORMATION			
Catalog No.	Description	Case Quantity	
303345	B-D® Blunt Plastic Cannula	100 per box	1,000 per case
303346	3cc Syringe with B-D® Blunt Plastic Cannula	100 per box	800 per case
303347	5cc Syringe with B-D® Blunt Plastic Cannula	100 per box	400 per case
303348	10cc Syringe with B-D® Blunt Plastic Cannula	100 per box	400 per case

**Interlink®**

MAKE THE CONNECTION

**BECTON  
DICKINSON**

Becton Dickinson Canada Inc. • 2771 Bristol Circle • Oakville, Ontario L6H 6R5  
B-D is a registered trademark of Becton Dickinson and Company.  
Interlink is a trademark of Baxter Healthcare Corporation.  
www.bd.com  
BDM 05/99-01