

Bovie®



ClearSpec® Speculum



* Patent US 6,432,048 B1

THE CHOICE IS CLEAR

- The only speculum clinically proven to provide complete visualization of the cervix in 3 times as many patients¹
- May help to improve the diagnosis by providing a more complete view of the cervix
- Sheath lateral ports facilitate tissue access for biopsy procedures
- Improves workflow while decreasing procedure time and patient discomfort due to less need for speculum repositioning

- Reusable LED Illuminator supplies two hours of light with automatic on/off switch to save battery life
- Illuminator conveniently fits into the handle of the speculum providing lighting without cords or recharging of batteries
- Illuminator is less expensive than traditional rechargeable light sources or illuminators that can be accidentally thrown away



* Patent Pending

ADVANCED SPECULUM DESIGN


Enabling improved visualization of the cervix.¹




ClearSpec®

PROBLEM

Standard competitors' specula only provide single-plane tissue retraction




View of the cervix using a standard competitors' speculum




SOLUTION

ClearSpec provides multi-planar tissue retraction



View of the cervix using ClearSpec



¹Hill et al., Sheathed versus standard speculum for visualization of the cervix. IJGO. 2014; 125:116-120

In this study the investigators were able to visualize the entire cervix in 42 (61.8%) patients when using the sheathed speculum [ClearSpec] compared with 11 (16.4%) patients undergoing standard speculum examination (p<0.0001).

The ClearSpec Products

MFG #	Item Description	Packaging Configuration
102	ClearSpec Single Use Vaginal Speculum with Illuminator - Small	25/bx
103	ClearSpec Single Use Vaginal Speculum with Illuminator - Medium	25/bx
104	ClearSpec Single Use Vaginal Speculum with Illuminator - Large	25/bx
105	ClearSpec LED Speculum Illuminator	12/bx

For more information, visit BovieMedical.com



Bovie Medical • Clearwater, FL
www.BovieMedical.com • Sales@BovieMed.com
 U.S. Phone 1 800 537 2790 • Int'l. Phone +1 727 384 2323

©2016 Bovie Medical Corporation, All Rights Reserved.
 Bovie® is a registered trademark of Bovie Medical Corporation.
 ClearSpec® is a registered trademark of ClearSpec LLC