

AARON<sup>®</sup> 950<sup>™</sup> -G

*Bovie*<sup>®</sup>

# ELECTROSURGERY FOR THE OFFICE AND SURGERY CENTER



## COMPLETE

The Aaron 950-G is packaged with the essential accessories needed to perform an OB/GYN LEEP procedure

## CONVENIENCE

The full system is situated on a convenient mobile stand for effortless portability

## AFFORDABLE

The most affordable OB/GYN system on the market

## PERFECTION




A variety of LLETZ loop electrodes are available for perfect surgical selection

## RELIABILITY

The Aaron 950-G is backed by Bovie's two (2) year manufacturer's limited warranty

ELECTROSURGICAL GENERATOR

AARON<sup>®</sup> 950<sup>™</sup>-G  
INCLUDES

3-Button Autoclavable Handpiece	Non-Sterile Sheath	Blunt & Sharp Dermal Tips
		
A901	A910 (10 Included)	A806 & A804 (50 each)

Disposable Electrodes	Solid Disposable Grounding Pad	Reusable Grounding Cord	Wall Mount Kit
			
ES01 (3 Included), ES02, ES20	A1202 (5 Included)	A1202C	A837

Smoke Shark <sup>™</sup> Smoke Evacuator & Pneumatic Footswitch	18-Hour Filter	6' Non-sterile Tube & Reducer Fitting	Mobile Stand & Tray
			
SE01	SF18	786T & SERF	A812-C

ESLK - Disposable Electrosurgical Loop Electrode Kit	ES07- Extended Ball Disposable Electrode - 5/box
	

Additional accessories available — contact your Bovie representative for details.

TECHNICAL  
SPECIFICATIONS

**A950 Technical Specifications...**

<b>Dimensions:</b>		<b>Line Voltage:</b>	<b>Line Frequency:</b>
Height	16.0 cm (6.3 in.)	A950: 120 VAC	50 - 60 Hz
Width	22.8 cm (9 in.)	A950-220: 240 VAC	
Depth	11.4 cm (4.5 in.)		
Weight	< 3.2 kg (< 7 lbs)		
<b>Monopolar</b>	<b>Output Power</b>	<b>Output Frequency</b>	<b>Rep Rate</b>
Cut	60 W @ 500 Ω	357 kHz ± 50 kHz	N/A
Blend	60 W @ 800 Ω	357 kHz ± 50 kHz	30 kHz ± 5 kHz
Coagulation	60 W @ 1000 Ω	375 kHz ± 50 kHz	60 kHz ± 5 kHz
Fulguration	35 W @ 1000 Ω	575 kHz ± 50 kHz	30 kHz ± 5 kHz
Bipolar	30 W @ 200 Ω	520 kHz ± 14 kHz + 50 kHz	19 kHz ± 5 kHz
<b>Warranty:</b> limited two (2) year on unit only			

**SE01 Technical Specifications...**

<b>Dimensions:</b> (without filter)	
Height	29.9 cm (11.8 in.)
Width	17.8 cm (7 in.)
Depth	46.4 cm (18.3 in.)
Weight	< 6.5 kg (< 14 lbs)
<b>Line Voltage:</b>	
SE01: 120 VAC	
SE01-220: 240 VAC	
<b>Line Frequency:</b>	
50 - 60 Hz	
<b>Warranty:</b> limited two (2) year on unit only	



Bovie Medical Corporation • 5115 Ulmerton Road • Clearwater, FL 33760 • www.boviemed.com • sales@boviemed.com  
U.S. Phone 1-800-537-2790 • Fax 1-800-323-1640 • Int'l Phone +1-727-384-2323 • Fax +1-727-347-9144

Bovie<sup>®</sup> and Aaron<sup>®</sup> are registered trademarks of Bovie Medical Corporation