

SMOKE SHARK™

Bovie®

SMOKE EVACUATION FOR THE OFFICE AND SURGERY CENTER



SMOKE EVACUATOR

EFFORTLESS

Simple setup and easy-to-operate

DISCREET

Quiet and effective smoke evacuation

AFFORDABLE

Long-life 18-hour four-stage filter keeps cost down per procedure and minimizes waste

EFFECTIVE

Excellent particle capture at the surgical site, down to 0.1 to 0.2 μm with 99.999% efficiency

RELIABILITY

The Smoke Shark is backed by Bovie's two (2) year manufacturer's limited warranty

STANDARD SETUP

Long-Life 18-hr Filter*	Pneumatic Footswitch	6' non-sterile tube*	Arm Stand*
			
SF18	FSSEP (included)	786T	SEAS

*sold separately

Customize your Smoke Shark™ for a variety of applications by choosing from the accessories below...

6' tube (sterile)	6' tube with wand & tip (non sterile)	6' tube with wand & tip (sterile)	Reducer fitting	Reducer fitting, sterile
				
786TS	SETW	SETWS	SERF	SERFS

Laparoscopic tubing kit (sterile)	Laparoscopic tubing kit with valve (sterile)	6' tubing for ValleyLab® & IEC Pencil (sterile)	Laser-resistant wand (non-sterile)	Laser-resistant wand (sterile)
				
SELKS	SELVS	SEVL	SELW	SELWS

Additional Accessories

Wall plate to hold canister	Fluid collection canister	Laser mask	Canister ring holder	remote switch activator
				
SEWP	SEFC	SELM	SERH	SERS

Technical Specifications...

Dimensions:		Line Voltage:	Line Frequency:
Height	29.9 cm (11.8 in.)	SE01: 120 VAC	50 - 60 Hz
Width	17.8 cm (7 in.)	SE01-220: 240 VAC	
Depth	46.4 cm (18.3 in.)		
Weight	< 6.5 kg (< 14 lbs) w/o filter		
Warranty: two (2) year			



OTHER OPTIONAL ACCESSORIES

TECHNICAL SPECIFICATIONS

