

Clinton Parallel Bars are made with the patient and therapist in mind. Sturdy, steel upright construction and smooth stainless steel handrails provide the strength for daily use. Simple height adjustment makes it easy on therapists who need to accommodate multiple individuals and exercise regimens.

Clinton Parallel Bars feature:

- One piece, 1 1/2" stainless steel, width adjustable handrails with end bumpers (3.8 cm)
- Uprights are ultra-heavy, square, chrome-plated, telescoping, steel tubing
- Powder-coated, all-steel, base tube upright
- Height is determined by recessed number strip
- All models, except wall mounted folding unit, feature width adjustment
- Platform models feature 1 1/2" thick, hardwood, plywood bases with ultra-durable, UV cured, satin finish (3.8 cm)
- Safety treads at both ends of platform models make it safer for patient to turn and exit wheelchair
- All models are for physical therapy use only

Platform Mounted Parallel Bars

- Stainless steel adjustable handrails
- Optional 045 Obduction Board is available (see page 41)
- 400 lbs. load capacity under normal use (181.43 kg)

3-2007	length	width range	height range
	7'	18" – 28"	26" – 39"
	2.13 m	45.72 – 71.12 cm	66 – 99 cm

3-2010 10' long model (3 m)

3-2012 12' long model (3.65 m)



3-2007