


RECOVERY COUCHES

Clinton Recovery Couches come in nearly as many styles and configurations as there are uses for comfortable, cost-effective couches. They're perfect for school nurse offices, outpatient surgery centers or anywhere temporary rest is needed for children or adults.

Recovery Couches feature:

- 7 versatile designs
- 2" firm foam padding (5 cm)
- Solid plywood, frame construction
- Easy clean laminate models, classic chrome or wood leg models
- Paper dispenser included with all models
- Premium, stain-resistant, woven, knit-backed vinyl upholstery
- 250 lbs. load capacity under normal use (113 kg)



3620

Hardwood Leg Couch

- Solid plywood, frame construction for strength and durability
- Non-adjustable, wedge headrest
- Solid, natural finish, hardwood legs
- Paper dispenser included
- 250 lbs. load capacity under normal use (113 kg)

3620-24	length	width	height
	72"	24"	18"
	183 cm	61 cm	46 cm

3620-27 27" width model (68.58 cm)



RECOVERY COUCHES

Chrome Leg Couch

- Solid plywood, frame construction
- Non-adjustable wedge headrest
- Chrome-plated, round, steel legs
- Paper dispenser included
- 250 lbs. load capacity under normal use (113 kg)

3600-24	length	width	height
	72"	24"	18"
	183 cm	61 cm	46 cm

3600-27 27" width model (68.58 cm)

Panel Leg Couch

- Steel frame construction
- Built-in, non-adjustable wedge headrest
- Easy-clean laminate legs
- Paper dispenser included (*mounted to head end leg*)
- 250 lbs. load capacity under normal use (113 kg)

3610-24	length	width	height
	72"	24"	18"
	183 cm	61 cm	46 cm

3610-27 27" width model (69 cm)

Upholstered Apron Couch with Double Drawer Storage

- Easy-clean, all laminate construction
- Adjustable, pillow wedge headrest
- Cabinet style base with flat panel laminate & steel drawers (available with optional 055 locks)
- Built-in paper dispenser
- 300 lbs. load capacity under normal use (136 kg)

3713-16	length	width	height
	72"	27"	20"
	183 cm	69 cm	51 cm

3713-15 with flat foam adjustable headrest

3713-10 with non-adjustable pillow wedge

Upholstered Apron Couch with Sliding Doors

- Easy-clean, all laminate construction
- Adjustable, pillow wedge headrest
- Cabinet style base with flat panel laminate & sliding doors (available with optional 055 locks)
- Built-in paper dispenser
- 300 lbs. load capacity under normal use (136 kg)

3770-16	length	width	height
	72"	27"	20"
	183 cm	69 cm	51 cm

3770-15 with flat foam adjustable headrest

3770-10 with non-adjustable pillow wedge



3600



3610



3713-16

Built-in paper dispenser



Optional Adjustable Wedge Headrest

016

- For models 3600, 3610 & 3620



Breathe easy with Clinton ECO PVC-Free upholstery.

Our soft, durable PVC-Free upholstery can be applied to ANY Clinton product that features traditional upholstery. Simply add the 117 optional suffix (*upcharge applies*) to your order number, and select any one of the eight earth hues from our ECO range of colors.

Plus, Clinton products utilize CARB II compliant plywoods and laminate panels, and our solid wood legs are top-coated with a water-based finish.

Clinton ECO makes sense for a wide range of reasons — for patients, for staff and everyone who cares about extra environmental safety.

TABLES/COUCHES



SEATING



PHLEBOTOMY

