

Wall Weight Bar Rac

- Plastic mounting hooks will not scratch bars
- Laminate surface
- Mounting holes included
- Holds 8 weight bars (Not included)

4716	width	height
	23 1/2"	36 1/2"
	59.69 cm	92.7 cm



4716

Floor Weight Bar Rac

- Easy-clean all laminate construction
- Holds 12 weight bars (Not included)

4717	width	depth	height
	15 7/8"	12"	13"
	40.32 cm	30.48 cm	33 cm



4717

4719

Mobile Weight Bar Rac

- Easy-clean all laminate construction
- 3" swivel casters (7.62 cm)
- Holds 12 weight bars (Not included)

4719	width	depth	height
	17"	12"	17 1/2"
	43.18 cm	30.48 cm	44.45 cm

Wall Weight Bar Rac

- Plastic mounting hooks will not scratch bars
- Laminate surface
- Mounting holes included
- Holds 5 weight bars (Not included)

4718	depth	height
	20"	5 1/2"
	50.8 cm	13.97 cm



Innovative grip system holds weight bars in place



4718

Wall Mount Ladder Bar Rac

- Can be used as a flexion shoulder ladder and/or a BarRac to store bars (Weight bars not included.)
- 8 levels provide a wide range of motion
- Sitting to standing exercise can be achieved by combining 2 wall mount units and mounting them vertically, one above the other.
- Virtually unbreakable polycarbonate rods

4716B	width	depth	height
	20 1/2"	36 1/2"	4 1/2"
	52 cm	90.17 cm	11.43 cm



4716B

Desktop Ladder Bar Rac

- Can be used as a flexion shoulder ladder and/or a BarRac to store bars (Weight bars not included.)
- 8 levels provide a wide range of motion
- Unit folds for easy storage
- Virtually unbreakable polycarbonate rods

4716D	width	depth	height
	20 1/2"	36 1/2"	4 1/2"
	52 cm	90.17 cm	11.43 cm

BarRacs 4716B & 4716D can also be used to store weight bars.