

Galaxy DualRac

- 5 independent, quick-change positions for holding your choice of exercise band (50 yds. (45.72 m), hollow core center only). All items shown not included.)
- Solid surface—no pegboard to chip
- Cuff weight “sliding hook system” makes lateral positioning easy
- Strong polycarbonate dumbbell rods are angled for added stability
- Protective bumper moulding around base ends
- 3" swivel casters/2 locking (7.62 cm)
- Holds up to 56 cuff weights, 24 dumbbells, 10 weight bars and 5 rolls of exercise band (Not included)

5119	width	depth	height
	36"	18"	65"
	91.44 cm	45.72 cm	165 cm



5119

Constellation DualRac

- 5 independent, quick-change positions for holding your choice of exercise band (50 yds. (45.72 m), hollow core center only). All items shown not included.
- Solid surface – no pegboard to chip
- Cuff weight “sliding hook system” makes lateral positioning easy
- Strong polycarbonate dumbbell rods are angled for added stability
- Protective bumper moulding around base ends
- 3" swivel casters/2 locking (7.62 cm)
- Holds up to 64 cuff weights, 16 dumbbells, and 5 rolls of exercise band (Not included)

5121	width	depth	height
	25 1/2"	18"	65"
	64.77 cm	45.72 cm	165 cm



5121

Back

Aurora DualRac

- Innovative design holds dumbbells and cuff weights
- Solid surface – no pegboard to chip
- Cuff weight “sliding hook system” makes lateral positioning easy
- Strong polycarbonate dumbbell rods are angled for added stability
- Protective bumper moulding around base ends
- 3" swivel casters/2 locking (7.62 cm)
- Holds up to 48 cuff weights and 22 dumbbells, or combination of dumbbells and weight bars (Not included)

5120	width	depth	height
	25 1/2"	18"	65"
	64.77 cm	45.72 cm	165 cm



5120

Back

Celestial DualRac with Mirror

- Mirror with ANSI 297 safety backing on reverse
- 5 independent, quick-change positions for holding your choice of exercise band (50 yds. (45.72 m), hollow core center only). All items shown not included.
- Solid surface – no pegboard to chip
- Cuff weight “sliding hook system” makes lateral positioning easy
- Strong polycarbonate dumbbell rods are angled for added stability
- Protective bumper moulding around base ends
- 3" swivel casters/2 locking (7.62 cm)
- Holds up to 16 cuff weights, 16 dumbbells, and 5 rolls of exercise band (Not included)

5121M	width	depth	height
	25 1/2"	18"	65"
	64.77 cm	45.72 cm	165 cm



5121M

Back

LAMINATE COLOR CHOICES



LAMINATE COLOR CHOICES

