


**Zeus SuperRac**

- 6 independent, quick-change positions for holding your choice of exercise band (50 yd. (45.72 m), hollow core center only). All items shown not included.
- Rac design uses all 4 sides for greater storage capacity
- Minimum floor space requirement
- Raised edge on top provides extra storage space
- Cuff weight "sliding hook system" makes lateral positioning easy
- Strong polycarbonate dumbbell rods are angled for added stability
- 4 drawers and open shelves for additional storage needs
- Protective bumper moulding around base
- 4" swivel casters/2 locking (10.16 cm)
- Holds up to 22 dumbbells, 20 weight bars, 4 (25 yds.) rolls of Velcro,® 6 rolls of exercise band and 48 cuff weights

5140	width	depth	height
	28"	28"	66"
	71.28 cm	71.28 cm	167.64 cm



**Atlas TowerRac with Mirror**

- 6 independent, quick-change positions for holding your choice of exercise band (50 yds. (45.72 m), hollow core center only). All items shown not included.
- Rac design uses all 4 sides for greater storage capacity
- Minimum floor space requirement
- Mirror has polished edges and ANSI 297 safety backing
- Raised edge on top provides extra storage space
- Cuff weight "sliding hook system" makes lateral positioning easy
- Strong polycarbonate dumbbell rods are angled for added stability
- Protective bumper moulding around base
- 4" swivel casters/2 locking (10.16 cm)
- Holds up to 48 cuff weights, 22 dumbbells, 20 weight bars and 6 rolls of exercise band

5139	width	depth	height
	28"	28"	66"
	71.28 cm	71.28 cm	167.64 cm



LAMINATE COLOR CHOICES



**Fraser SlimRac**

- 5 independent, quick-change positions for holding your choice of exercise band (50 yds. (45.72 m), hollow core center only). All items shown not included.
- Rac design uses all 4 sides for greater storage capacity
- Minimum floor space requirement
- Cuff weight "sliding hook system" makes lateral positioning easy
- Strong polycarbonate dumbbell rods are angled for added stability
- Protective bumper moulding around base
- 4" swivel casters/2 locking (10.16 cm)
- Holds up to 64 cuff weights, 27 dumbbells, 12 weight bars and 5 rolls of exercise band

5133A	width	depth	height
	28"	22"	65 1/2"
	71.28 cm	55.88 cm	166.37 cm



**Corbin SlimRac**

- 5 independent, quick-change positions for holding your choice of exercise band (50 yds. (45.72 m), hollow core center only). All items shown not included.
- Rac design uses all 4 sides for greater storage capacity
- Minimum floor space requirement
- Cuff weight "sliding hook system" makes lateral positioning easy
- Strong polycarbonate dumbbell rods are angled for added stability
- Protective bumper moulding around base
- 4" swivel casters/2 locking (10.16 cm)
- Holds up to 88 cuff weights, 27 dumbbells and 5 rolls of exercise band

5133	width	depth	height
	28"	22"	65 1/2"
	71.28 cm	55.88 cm	166.37 cm



LAMINATE COLOR CHOICES

