

**Mobile Adult Mirror**

- 1/4" thick distortion-free glass (0.63 cm) mirror with ANSI 297 safety backing
- Hardwood, plywood frame
- Rolls easily on 3" swivel casters (7.62 cm)
- Finished laminated back
- Black tee-moulding bumper strips on base sides

6210	width	depth	height
	26 1/2"	21"	72"
	67.31 cm	53.34 cm	182.88 cm

**Mobile Mirror with Cuff Weights**

- Same as model 6210 with cuff weights on reverse side
- Holds up to 48 cuff weights (Not included)

6210-9	width	depth	height
	26 1/2"	21"	72"
	67.31 cm	53.34 cm	182.88 cm

**Mobile Youth Mirror**

- Smaller version of model 6210 (see 6210 photo)

6222	width	depth	height
	17 1/2"	21"	57"
	44.45 cm	53.34 cm	144.78 cm



**Wall Mounted Mirror**

- 1/4" thick distortion-free glass mirror with ANSI 297 safety backing (0.63 cm)
- Hardwood, plywood frame
- Easily attaches to wall with top and bottom mounting cleats

6220	width	depth
	25"	67"
	63.5 cm	170.18 cm

**Mobile Tri-Panel Mirror**

- Three 1/4" thick distortion-free glass mirror with ANSI 297 safety backing (0.63 cm)
- Hardwood, plywood frame
- Rolls easily on 3" swivel casters (7.62 cm)
- Finished laminated back
- View angles easily adjust

6223	open width	panel width	height
	60"	25"	70"
	152.4 cm	132.08 cm	177.8 cm

