

Clinton designs and manufactures pediatric/infant blood drawing stations that are easy to use, durable and patient friendly. Choose from sleek classic style models or Clinton Kids™ Imagination Series models with fun graphics embedded right in the laminate so they won't tear or peel off.

**Pediatric/Infant Series, Blood Drawing Stations feature:**

- Easy-clean padded tops
- Stainless Steel frame or durable laminate panel construction
- Nylon glides on all models
- 200 lbs. load capacity under normal use (90.71 kg)

**IMPORTANT NOTE:** Clinton Blood Drawing Stations are not engineered to restrain a fainting or falling occupant. Children and infants should NEVER be left unattended on or around the station.



6077



67114



67110

**Stainless Steel, Infant Blood Drawing Station**

- All-welded, stainless steel frame
- Top has 2 rims
- 2 safety straps
- Plastic work tray Included
- Removable middle padded section

**6077**

overall	width	length	height
	20 1/2" 52 cm	39" 99 cm	39 1/2" 100.33 cm

top	width	length
	20" 50.8 cm	32" 81.28 cm

**Panel Leg, Infant Blood Drawing Station with Shelf**

- All laminate construction
- Top rimmed on 3 sides for safety
- Easy-clean, center pad

**67114**

overall	width	length	height
	25" 63.5 cm	43" 109.2 cm	39 1/2" 100.33 cm

top	width	length
	25" 63.5 cm	43" 109.2 cm

**Panel Leg, Infant Blood Drawing Station**

- All laminate construction
- Top rimmed on 3 sides for safety
- Easy-clean, center pad

**67110**

overall	width	length	height
	25" 63.5 cm	43" 109.2 cm	39 1/2" 100.33 cm

top	width	length
	25" 63.5 cm	43" 109.2 cm

**Imagination Series/Rainforest Follies, Infant Blood Drawing Station with 2 Doors**

- All laminate construction
- Top rimmed on 3 sides for safety
- Easy-clean, center pad
- Large, 2-door storage
- Clinton Kids™ graphics embedded in the laminate will not tear or peel

**67132**

overall	width	length	height
	25" 63.5 cm	43" 109.2 cm	39 1/2" 100.33 cm

top	width	length
	25" 63.5 cm	43" 109.2 cm



67132

**Imagination Series/Space Place, Infant Blood Drawing Station with 2 Doors**

- All laminate construction
- Top rimmed on 3 sides for safety
- Easy-clean, center pad
- Large, 2-door storage
- Clinton Kids™ graphics embedded in the laminate will not tear or peel

**67135**

overall	width	length	height
	25" 63.5 cm	43" 109.2 cm	39 1/2" 100.33 cm

top	width	length
	25" 63.5 cm	43" 109.2 cm



67135

**Imagination Series/Ocean Commotion, Infant Blood Drawing Station with 2 Doors**

- All laminate construction
- Top rimmed on 3 sides for safety
- Easy-clean, center pad
- Large, 2-door storage
- Clinton Kids™ graphics embedded in the laminate will not tear or peel

**67136**

overall	width	length	height
	25" 63.5 cm	43" 109.2 cm	39 1/2" 100.33 cm

top	width	length
	25" 50.8 cm	43" 109.2 cm



67136