

Laminate Finger/Shoulder Ladder

- 36 smooth round corner steps
- Color coded every sixth step
- Easy to clean and sanitize gray surface
- Pre-drilled mounting holes

4758A	length	width
	47"	6"
	119.37 cm	15.23 cm

Wood Finger/Shoulder Ladder

- 36 vertical finger steps
- Solid wood with heavy topcoats
- Pre-drilled mounting holes

7220	length	width
	58"	2 1/2"
	147.51 cm	6.35 cm

Shoulder Wheel

- 2 chrome-plated height adjustment rails
- Handle adjusts to increase diameter range
- Adjustable resistance knobs
- 2 knobs adjust height in 1 1/2" increments (3.8 cm)

7213	depth	width	diameter range
	8"	38"	10" – 39"
	20.31 cm	96.52 cm	43.18 cm – 99 cm

Exercise Skate

- Padded Dove Gray vinyl top
- Top is mounted on 4 casters for easy movement
- Velcro® strap(s) secures hand


7206	length	width	height
Small Skate	5 1/2"	5 1/2"	3"
	13.97 cm	13.97 cm	7.62 cm

7212	length	width	height
Large Skate	12"	6"	3"
	30.47 cm	15.24 cm	7.62 cm

Powder Board Table

- Oval shape with patient cutout
- White high pressure laminate top with tee-moulding edges
- 2 knobs adjust height in 1 1/2" increments (3.8 cm)
- Uprights are ultra-heavy, square, telescoping, steel tubing with self-indexing, spring plunger (no holes to line up)
- Second knob removes any side-to-side play
- Height is determined by numbered height indicator strips (no holes to count)

75-23K	length	width	height range
	63"	40"	23" – 34"
	160 cm	101.6 cm	58.42 cm – 86.36 cm



4758A 7220 7213



7206 & 7212



75-23K

Nesting Stools

- All hardwood plywood construction
- Safety tread on all tops
- 350 load capacity under normal use

7320	length	width	height
set 1	12"	10"	2"
	30.48 cm	25.4 cm	5.08 cm
set 2	14"	12"	4"
	35.56 cm	30.48 cm	10.16 cm
set 3	18"	14"	6"
	45.72 cm	35.56 cm	15.24 cm
set 4	20"	16"	8"
	50.8 cm	40.64 cm	20.32 cm



7320

Adjustable Balance Platform

- Adjustable height platform mounted handrails
- Platform carpeted for safety
- Easy-grip 1 1/2" diameter stainless steel rails (3.8 cm)
- 350 load capacity under normal use

7360	length	depth	height range
	36"	36"	26" – 39"
	91.44 cm	91.44 cm	66.04 cm – 99.06 cm



7360

Ramp and Curb Training Set

- Adjustable unit/remove and sequence curbs
- Curbs store beneath platform
- All hardwood plywood construction
- Safety tread on all stairs
- Reinforced with steel angle
- 350 load capacity under normal use

7330	length	width	height
ramp	84"	44"	8"
	213.36 cm	111.76 cm	20.32 cm
curb	32"	32"	2"
	81.28 cm	81.28 cm	5.07 cm
curb	36"	36"	4"
	91.44 cm	91.44 cm	1.16 cm
curb	40"	40"	6"
	101.6 cm	101.6 cm	15.24 cm



7330

Quad Board

- 10 position, adjustable incline
- Gray, padded surface folds flat for storage
- Natural finish, wood construction
- 350 load capacity under normal use

7370	length	depth	height range
	29"	5 1/2"	7" – 13"
	73.66 cm	13.97 cm	17.78 cm – 33 cm



7370

Floor Style Powder Board Table

- Smooth, light colored top
- 4 stationary legs and full frame
- Folds flat for easy storage

7350	length	width	height
	40"	29"	7"
	101.6 cm	73.66 cm	17.78 cm



7350