

## TREATMENT CABINETS

### Laminate Top, Mobile, Equipment Cabinet

- Mounted on a stable 5-leg black nylon base
- 2 self-closing, Euro-style drawers with bumpers
- Laminate top rimmed on 3 sides
- All laminate cabinet and top

8920	width	depth	height
	22"	17"	33"
	55.88 cm	43.18 cm	83.82 cm

### Laminate Top, Mobile, Equipment Cart with 1 Drawer and 1 Shelf

- Mounted on a stable 5-leg black nylon base
- 1 Euro-style drawer
- Top and back have wire feed cut-outs
- Top rimmed on 3 sides
- All laminate cabinet

8940	width	depth	height
	22"	17"	33"
	55.88 cm	43.18 cm	83.82 cm

### Moulded Top, Mobile, Equipment Cabinet

- One-piece, vacuum formed gray haircell textured top rimmed on all sides
- Mounted on a stable 5-leg black nylon base
- 2 self-closing, Euro-style drawers with bumpers
- Top rimmed on 4 sides
- All laminate cabinet
- Ships knocked-down

8923-A	width	depth	height
	24 <sup>1</sup> / <sub>2</sub> "	18 <sup>3</sup> / <sub>8</sub> "	33"
	62.22 cm	46.67 cm	83.82 cm

### Optional Basket

- Wire basket with white epoxy finish

089	width	depth	height
	17 <sup>1</sup> / <sub>2</sub> "	6 <sup>3</sup> / <sub>4</sub> "	5"
	44.45 cm	17.14 cm	12.7 cm



8920



8940



8923-A



089

## HAND THERAPY TABLES

### Mobile Hand Therapy Table

- Spring loaded height adjustment with 2 hand control levers
- 1<sup>1</sup>/<sub>8</sub>" Gray laminate top only (2.85 cm)
- Chrome base
- Removable casters for stationary application

74-10S	length	width	height range
	32"	18"	28"- 43"
	81.28 cm	45.72 cm	71.12 cm-109.22 cm

### Hand Therapy Table\*

- Height adjustable legs with self indexing spring loaded plunger
- Easy read, numbered height strip

74-15K*	length	width	height range
	36"	19"	23"- 33"
	91.44 cm	42.26 cm	58.42 cm-83.82 cm

### Soft Curve Gas Lift Table\*

- 1<sup>1</sup>/<sub>8</sub>" laminate top with sit-close cut-out on all sides (2.85 cm)
- Gas spring assist lever for easy height adjustment
- 30" (76.2 cm) base also available with optional casters (-081)

74-17G*	length	width	height range
	36"	36"	27"- 38"
	91.44 cm	91.44 cm	68.58 cm-96.52 cm

### Soft Curve Hydraulic Lift Table\*

- 1<sup>1</sup>/<sub>8</sub>" laminate top with sit-close cut-out on all sides (2.85 cm)
- Foot pump height adjustable base with swivel snug-top
- 28" chrome plated steel base (20.32 cm)

74-14H*	Length	width	height range
	36"	36"	28"- 36"
	91.44 cm	91.44 cm	71.12 cm-91.44 cm

### Hand Therapy Station with Shelf Unit

- 1<sup>1</sup>/<sub>8</sub>" laminate top (2.85 cm)
- 3 shelf storage
- Natural finish only

74-110	length	depth	height
overall	45"	20"	31"
	114.3 cm	50.8 cm	78.74 cm
shelf	20"	12"	30 <sup>1</sup> / <sub>2</sub> "
	50.8 cm	30.48 cm	77.47 cm

74-120 Same as 74-110 without shelf unit



74-10S



74-15K



74-17G



74-14H



74-110

74-120 Station without shelf unit

\* Specify Maple, Gray or Natural top colors when ordering