



“Green” isn't just a catch phrase. It's a choice. At Clinton, we take environmentally-friendly products seriously by designing them to provide additional safeguards for indoor air quality. We use materials containing no added urea-formaldehyde, use recyclable materials wherever possible, and use PVC-free vinyl on all our Clinton-Eco products. Overall, Clinton provides a broad range of eco-friendly equipment that looks great and meets your “green” expectations.

Clinton-Eco, Wood Tables feature:

- Clinton uses CARB I and CARB II compliant plywood and laminate panels with no added urea-formaldehyde
- Water-based or UV-cured topcoats
- PVC-free upholstery
- PVC-free, edge tape for shelving

Clinton-Eco, Metal Tables feature:

- 100% recyclable metal table frame
- Eco-friendly applied VOC-free, powder-coated metal
- PVC-free upholstery
- PVC-free edge tape for shelving
- Clinton uses CARB I & CARB II compliant, laminated panels with no added urea-formaldehyde

Clinton-Eco, Cabinet Tables feature:


- Clinton uses CARB I and CARB II compliant plywoods and laminated panels with no added urea-formaldehyde
- PVC-free upholstery
- PVC-free, edge tape for shelving
- 100% recyclable metal components



CLINTON-ECO PVC-FREE UPHOLSTERY



CLINTON-ECO BASE FINISHES



Due to the color limitations of the printing process, actual colors may vary. Actual laminate and upholstery samples are available by request.



Eco, Wood, Treatment Table

- Features all hardwood legs, H-bracing and adjustable backrest
- Rounded corner top
- 400 lbs. load capacity under normal use (181.43 kg)
- Also available with optional flat top (088)*

81010-24	length	width	height
	72"	24"	31"
	182.88 cm	60.95 cm	78.73 cm

81010-27	27" width model (68.58 cm)		
----------	----------------------------	--	--

81010-30	30" width model (76.2 cm)		
----------	---------------------------	--	--



81010

Eco, Wood, Treatment Table with Drawers

- Features all hardwood legs, laminate shelf, two drawers and adjustable backrest
- Rounded corner top
- 400 lbs. load capacity under normal use (181.43 kg)
- Also available with optional flat top (088)*

81013-27	length	width	height
	72"	27"	31"
	182.88 cm	68.57 cm	78.73 cm

81013-30	30" width model (76.2 cm)		
----------	---------------------------	--	--



81013

Eco, Wood, Treatment Table with Shelf

- Features all hardwood legs, full laminate shelf and adjustable backrest
- Rounded corner top
- 400 lbs. load capacity under normal use (181.43 kg)
- Also available with optional flat top (088)*

81020-24	length	width	height
	72"	24"	31"
	182.88 cm	60.95 cm	78.73 cm

81020-27	27" width model (68.58 cm)		
----------	----------------------------	--	--

81020-30	30" width model (76.2 cm)		
----------	---------------------------	--	--



81020

Eco, Wood, Treatment Table with Shelving

- Features all hardwood legs, laminate shelving and adjustable backrest
- Rounded corner top
- 400 lbs. load capacity under normal use (181.43 kg)
- Also available with optional flat top (088)*

81030-24	length	width	height
	72"	24"	31"
	182.88 cm	60.95 cm	78.73 cm

81030-27	27" width model (68.58 cm)		
----------	----------------------------	--	--

81030-30	30" width model (76.2 cm)		
----------	---------------------------	--	--



81030

* See all options on page 25