

Ball Base

- Provides stability for beginners
- Helps keep ball from rolling
- Cut-outs in base allow for band exercise while seated on the ball
- Also can be used to store balls

8402	diameter
	45 cm

Squeeze Balls

- Excellent for hand exercise
- 4 different densities for progressive exercise
- Non-stick material
- 50 mm diameter

	density	color
8301	soft	orange
8302	medium	green
8303	firm	teal
8304	extra firm	purple

Ball Stackers

- Flexible adjustable plastic ring allows the stacking of 3 exercise balls
- Sold as a set of 3 rings

8405	accommodates
	45 cm – 85 cm balls

Double Level Ball Storage Rack plus Balls

- All Gray laminate construction
- Store up to 6 therapy balls
- Easy assembly with steel fasteners
- Ships KD by FedEx® Ground or UPS®
- Set of heavy duty, burst-resistant exercise balls
1 each of—30 cm • 45 cm • 55 cm • 65 cm • 75 cm • 85 cm • plus air pump & inflation tape

8088	length	depth	height
	62"	19 1/2"	30"
	157.48 cm	49.53 cm	76.2 cm


8008	rack only
-------------	-----------

Triple Level Ball Storage Rack plus Balls

- All Gray laminate construction
- Store up to 9 therapy balls
- Easy assembly with steel fasteners
- Ships KD by FedEx® Ground or UPS®
- Set of heavy duty, burst-resistant exercise balls
2 each of—30 cm • 45 cm • 55 cm
1 each of— 75 cm • 85 cm • plus air pump & inflation tape

8099	length	depth	height
	62"	19 1/2"	60"
	157.48 cm	49.53 cm	152.4 cm

8009	rack only
-------------	-----------



NOTE: Colors of individual balls are subject to change.

Heavy-Duty Exercise Balls*

- Bright colored, heavy-walled, burst-resistant balls
- Deflates slowly if punctured*
- Made from durable, ribbed PVC material
- Shipped in polybags
- For indoor use on a clean, dry surface only

* Deflation rate is determined by size of puncture and weight of person. A ball inflation tape is included to achieve and maintain correct ball diameters.

	color	inches	centimeters	patient height
8030	green	15 1/2"	39 cm	–
8045	yellow	17 1/2"	44 cm	55" – 60"
8055	red	21 1/2"	55 cm	61" – 66"
8065	blue	25 1/2"	65 cm	67" – 71"
8075	green	29 1/2"	75 cm	72" – 75"
8085	black	33 1/2"	85 cm	76" – 80"



Heavy-Duty Exercise Balls

Peanut Balls*


- Burst resistant, latex-free balls
- Deflates slowly if punctured*
- Made from durable, ribbed, heavy-walled, PVC material
- Peanut shape helps control motion in the center of ball
- For inside use only

* Deflation rate is determined by size of puncture and weight of person.

	size	color
8540	40 cm	yellow
8550	50 cm	red
8560	60 cm	blue
8570	70 cm	green



Peanut Balls



7801

Thin Exercise Mat

- Designed to help stabilize feet and hands during exercise

NOTE: Colors of individual exercise mats are subject to change.

7801	length	width	thickness
	68"	36 1/2"	1/8"
	172.72 cm	92.71 cm	.317 cm

Ball Inflation Tape**

- Used to measure the outside diameter of the exercise balls during the inflation process. The tape easily attaches to the inflation needle or the plug. This is highly recommended to make measuring the diameter of the balls easy during inflation.
- Included with all exercise balls, but also sold separately



8406



8403

8406	
-------------	--

8403	Air Pump to Inflate Balls
-------------	---------------------------

*NOTE: Colors of individual balls are subject to change.