

INTEGRAL STEP STOOL TABLES

NOTE: These products will be available for order after September 1, 2016

Clinton, integral step stool, treatment tables are truly a step up from regular treatment tables. With an easy access, pull-out, welded steel, step stool built right into the table, these Clinton tables put convenience right where it's needed. Plus, the stool stores neatly under the table when not in use, saving valuable exam room space.

All Integral Step Stool Tables feature:

- Built-in, pull-out, welded steel, step stool
167/8" L x 191/2" W x 73/4" H
(42.84 cm L x 49.53 cm W x 19.68 cm H)
- Supply storage features
- Paper dispenser standard
- Available in woodtones and solid color laminates
- Wide selection of upholstery colors
- 400 lbs load/force capacity under normal use



Integral Step Stool, Treatment Table

- All laminate, panel leg design
- Handy, welded steel, pull-out, step stool facing front at the foot end of the table
- 2-drawer storage
- 6-position, adjustable backrest
- 400 lbs. load capacity under normal use (181.43 kg)

9402	length	width	height
	72"	27"	32"
	193 cm	68.58 cm	81.28 cm



9402



8801

Integral Step Stool, Treatment Table

- All laminate, cabinet design
- Handy, welded steel, pull-out, step stool at the foot end of the table
- Adjustable, pneumatic backrest
- 2-door and 2-drawer convenient storage
- Doors can be on either side of table (specify when ordering)
- 400 lbs. load capacity under normal use (181.43 kg)

8801	length	width	height
	72"	27"	32"
	193 cm	68.58 cm	81.28 cm

Integral Step Stool, Treatment Table

- Easy-clean, all laminate construction
- Handy, welded steel, pull-out, step stool facing front at the foot end of the table
- Handy 1-door and 2-drawer storage
- Adjustable, pneumatic backrest
- 400 lbs. load capacity under normal use (181.43 kg)

9407	length	width	height
	72"	27"	32"
	193 cm	68.58 cm	81.28 cm



9407



9405

Integral Step Stool, Treatment Table

- All laminate, cabinet design
- Handy, welded steel, pull-out, step stool facing front at the foot end of the table
- Adjustable, pneumatic backrest
- 2-door and 2-drawer convenient storage
- 400 lbs. load capacity under normal use (181.43 kg)

9405	length	width	height
	76"	27"	32"
	193 cm	68.58 cm	81.28 cm

NOTE: This product will be available for order after September 1, 2016



Stool in Stored Position



Pneumatic backrest features extra table paper storage compartments