

Endo Dissect™

5 mm laparoscopic hand instrument



PRODUCT INFORMATION

Code	Description	Qty Per Case
174505	Endo Dissect™ Short 5 mm, Single Use Instrument with Monopolar Cautery	1
176645	Endo Dissect™ 5 mm, Single Use Instrument with Monopolar Cautery	6

MATERIALS OF MANUFACTURE

Code	
174505	301, 302, 303, 304, and 416 Stainless Steel, 10% Glass Filled Polycarbonate, Polycarbonate 20% glass, PVC GRADE 6K10, STYLE VW-1, Aluminum 6061-T6, Union Carbide
176645	301, 302, 303, 304, and 17-4 P.H. Stainless Steel, Polycarbonate 20% glass, 10% Glass filled polycarbonate, Aluminum 2011-T3, PVC Grade 6K10

SHIPPING INFORMATION

	Inner Box Size (mm)	Shipper Box Size (mm)
174505	L 162 × W 29 × H 448	L 184 × W 165 × H 451
176645	L 610 × W 191 × H 70	N/A

ADDITIONAL INFORMATION

Latex Free	✓
DEHP Free	✓
Sterilization method	Radiation
Shelf Life	5 years
CE Marked	Class IIb Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 17-weu-endo-dissect-tech-sheet-1516046

Medtronic
Further, Together