



COVIDIEN

positive results for life™



Multi-Purpose Sharps Container

Multi-Purpose Sharps Containers

Order Information

Order Code	Description	Ship Case
8920SA	.5 Gallon Needle Syringe Collection Container 4.5" H x 4.75" D x 4.75" W	60
8950SA	5 Quart Container, Round, Red 6.75" H x 8.75" D	40
8970	2 Gallon Red Container with Rotor Lid 10" H x 7.25" D x 10.5" W	20
8970AC	2 Gallon Red Container with Rotor Lid 10" H x 7.25" D x 10.5" W	5
8979C	2 Gallon Clear Container with Rotor Lid 10" H x 7.25" D x 10.5" W	20
8990SA	2 Gallon Red Container with Clear, Hinged Lid 10" H x 7.25" D x 10.5" W	40
8996	2 Gallon Red BIOMAX Container with Hinged, Sealing Gasket Lid 10"H x 7.25"D x 10.5"W	20
8535SA	2 Gallon Clear Container with Rotor/Hinged Lid 10" H x 7.25" D x 10.5" W	20
8969	2 Gallon Beige Container with Clear Lid 10" H x 7.25" D x 10.5" W	20
8527	3 Gallon Clear Container with Hinged, Rotor Lid 13.75"H x 6"D x 13.75"W	10
8527R	3 Gallon Red Container with Hinged, Rotor Lid 13.75"H x 6"D x 13.75"W	10
8527RAC	3 Gallon Red Container with Hinged, Rotor Lid 13.75"H x 6"D x 13.75"W	5

Features and Benefits

Features	Benefits
Containers lock for final disposal	To secure contents
Nestable	Saves valuable storage space
Adjustable rotor opening or hinged lids	Accomodate various sizes of sharps

Frequently Asked Questions

Question	Answer
Are the containers autoclavable?	The containers without gasket are autoclavable.

COVIDIEN, COVIDIEN with Logo, and TM marked brands are trademarks of Covidien AG or an affiliate.
^{TM*} - Third party trademark.

© 2012 Covidien. All rights reserved.