

PROMED®



Models 6127 / 6129

DIGITAL HEALTHCARE SCALES

- Large 500lb/225kg Capacity
- USA Made
- Body Mass Index




Detecto's ProMed® 6129 medical scale with height measurement rod is the right choice for home, health, and clinical office weighing. This all-USA-manufactured scale is designed for years of dependable and accurate weighing, plus it's versatile enough to meet all clinical health requirements. With a capacity of up to 500 lbs, the 6129's scale platform is designed to monitor weight in the home or medical office. Prescribed by dietitians and physicians, the inclusive BMI (Body Mass Index) feature is easy to read for the monitoring of personal health and medical needs.

*Detecto
Model 6129
Digital
Healthcare Scale*

WAIST-HIGH DIGITAL HEALTHCARE SCALES

The all-steel construction of Detecto's ProMed® scales ensures years of reliable weighing and the baked-on epoxy paint finish is easy to clean for maintaining a sanitary environment in the home or office. Additionally, the large weighing platform adds stability and comfort for everyday use. Detecto is the original producer of the physician scale most seen in the medical industry and featured on the popular "How It's Made" television show.



Model 6127

MODEL	CAPACITY	HEIGHT ROD	SHIPPING WEIGHT
6127	500 lb x 0.2 lb / 225 kg x 0.1 kg	none	31 lb / 14 kg
6127KG	225 kg x 0.1 kg	none	31 lb / 14 kg
6129	500 lb x 0.2 lb / 225 kg x 0.1 kg	30" - 78" / 76 cm - 200 cm	34 lb / 15 kg
6129KG	225 kg x 0.1 kg	30" - 78" / 76 cm - 200 cm	34 lb / 15 kg

Platform: 10.75" W x 14.75" D / 27 cm W x 37.5 cm D
Display: Six-digit, seven-segment, 0.75" / 19 mm high LCD
Power: 6 "AA" size Alkaline, Ni-Cad or NiMH batteries (not included)
 OR an optional 100 to 240 VAC 50/60Hz 12 VDC 1A wall plug-in UL/CSA listed AC power supply (not included)



Model 6129KG

- Easy BMI calculation.
- Heavy-duty, USA-built, all-steel manufactured product for long-lasting durability.
- Large display for viewing of weight and BMI.
- Spacious weighing platform.
- Supports Alkaline or NiCad rechargeable batteries and optional AC adapter.
- 500 lb / 225 kg x 0.2 lb / 100 g, high-capacity digital scale.

OVERALL DIMENSIONS

Width: 10.875" / 28 cm
Height: 37.875" / 96 cm
Depth: 20.875" / 53 cm



Extra-large, 0.75 inch / 19 mm high digits on an easy-to-read LCD display.

Simply touch the BMI key and utilize the toggle arrows to adjust the patient height for BMI calculation.

Whether you require a scale for personal health or for the busy medical office, Detecto's ProMed® series scales are the models that best support any weighing needs you may have. ProMed® scales are fully supported by Detecto's complete family of customer service for questions, assistance and product support.

Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.



A Division of Cardinal Scale Manufacturing Co.

203 E. Daugherty, Webb City, MO 64870 USA
 Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001
 www.Detecto.com

Sold By:

