



**PRODUCT SPECIFICATIONS**



## EURO CHAIR SCALE

**MODEL: 6875  
PORTABLE  
ELECTRONIC**

Capacity: 400 lb x 0.2 lb/180kg x .1 kg  
Readout: Large easy-to-read LCD display  
Functions:

- Auto sleep mode - display goes off within two minutes after Lock/Weight key is depressed and the patient leaves the chair (lengthens battery life)
- Auto power off (after nine minutes of non-use)
- Lock/release of display
- Automatic/manual zero
- Motion detector

Key Operations: On/zero • Off • Weight/lock release • Body Mass Indicator Calculator

Color: White with black seat and

armrest  
Power: 6 "C" size Alkaline batteries or optional AC adapter.  
Features: Lightweight, easily moved on built-in wheels. Flip-up footrest for ease of patient access. Safety stops to prevent accidental tipping. Easily passes through a 30" door. Equipped with retractable armrests for ease in transferring patients to chair scale and to accommodate patients with casts.

Ship Weight: 70 lb



Cardinal Scale reserves the right to improve, modify, or enhance features and specifications without prior notice.

