



DIGITAL STADIOMETER

model **DHRWM**

Standalone Wall-Mount Digital Height Rod



Detecto's exclusive digital height technology. Receive height readings quickly, accurately, and comfortably.

Measure up to 6'7" / 200 cm

Automatic height viewing on the remote display, positioned where you need it with Detecto's wall-mount or desktop bracket.





ProDoc®

DIGITAL STADIOMETER model **DHRWM**

Standalone Wall-Mount Digital Height Rod

Features:

- Fast, accurate digital height LCD readouts
- Modern, clean design
- Display can be mounted on wall or desk
- View height in inches or centimeters
- Measures up to 6'7" / 200 cm

Specifications:

DISPLAY: Five-digit, seven-segment, 1 in (25 mm) high LCD

MEASURES: 43 in / 110 cm - 79 in / 200 cm

POWER: 6 "AA" Alkaline batteries (not included)

CABLE LENGTH: 69 in / 175 cm (*display to height rod*)

SHIPPING WEIGHT: 8.2 lb / 3.7 kg

Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.

Sold By:



DETECTO®
A Division of Cardinal Scale Manufacturing Co.
203 E. Daugherty, Webb City, MO 64870 USA
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001
www.Detecto.com