

Model DS4900

# Neonatal/Infant Scale Cart

Constructed of heavy gauge painted steel, the scale cart provides a convenient and safe way to move the neonatal or infant scale around a busy medical office or hospital pediatric ward.

The base of the cart features high impact plastic bumpers and a low center of gravity to prevent the cart from tipping over. Moving the scale is easy thanks to the 5" diameter wheels and ergonomic push handles. Two of the easy roll casters lock to ensure that the scale cart stays in place while the scale is being used.

The roomy shelves has plenty of space for medical supplies and the 3" deep drawer locks for security.

## CONSTRUCTION

ABS plastic and painted steel

## DIMENSIONS

Cart 29" W x 24.5" D x 37.5" H

Top Tray Dimension 25.75" L x 16.75" W x 0.5" D

## WARRANTY

**5** Five Years



The DS4900 accommodates either the DS4300 Neonatal Scale or the DS4100 Infant Scale.



**Doran Scales, Inc.** 1-800-264-4106  
1315 Paramount Parkway Fax: 630-879-0073  
Batavia, IL 60510 [www.doranmedical.com](http://www.doranmedical.com)