

Model DS6150

Remote Indicator Scale



The affordably priced Model DS6150 Remote Indicator Scale satisfies a variety of applications with a versatile design. Whether your application is a physician's office, acute care or visiting home nurse, the DS6150 can meet your needs. The indicator can be placed on a tabletop or attached to a wall using the mounting bracket. The DS6150 scale base is constructed of heavy gauge steel ensuring a long life and the complete scale system is backed by a three year warranty. The surface of the scale base is covered edge to edge with a non-skid material, ensuring patient safety.

Have the patient step onto the low-profile DS6150 platform and an accurate weight up to 500 pounds will be held on the display. The Weight Lock feature, a Doran exclusive, allows an accurate weight to be displayed regardless of a patient's movements. The Reweigh feature allows a second weight to be acquired while the patient remains on the platform.

Standard USB and RS-232 communications ports allow the recording of a patient's weight and BMI value by an EMR/EHR system.

MAXIMUM CAPACITY

500 x 0.2 lb (225 x 0.1 kg)

ACCURACY

+/- 0.6 lb (0.3 kg)

DIMENSIONS

Platform 15" W x 12" D x 1.5" H
Indicator 8" H x 6.5" W x 2" D

DISPLAYED UNITS

lb, kg

DISPLAY

1" High LCD

CONTROLS

Units, Reweigh, BMI, Up, Down, and Zero

BMI CALCULATOR

Input Data in FT/IN or CM

CONSTRUCTION

Indicator - Plastic
Platform - Steel

POWER

110 VAC Adapter - Included
(ULc Approved)
or Six AA Batteries

BATTERY LIFE

25,000 Weighments

WARRANTY

Three Years

OPTIONS

- RS-232 Communications Cable
- USB Communications Cable
- Welch Allyn® Spot Vital Signs® LXi Communication Cable
- AC Adapter
(110 - 220 VAC - CE Approved)
Available with Type A, C, G, or I plug for international use.

Model DS6150-WIFI



doranmedical.com
800-264-4107