

NEW!



Available
1 Pack

300 LB
WEIGHT CAPACITY

LiftWalker

LWWSW Lift Walker, Two Button with Assist Bars, 1/cs

FEATURES:

- For use by individuals who need assistance with sitting/standing on their own.
- Safer than being pulled by an attendant or caregiver.
- Reduces risk of injury to an attendant or caregiver.
- The retractable bars are easy to deploy.
- The Lift Walker makes it easy to stand because the individual uses biceps, chest, back, stomach and leg muscles all at once.
- Easy to lock/unlock foot brake.
- Easy push button mechanism may be operated by fingers, palms, or side of hand.
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing.
- Comfortable, soft foam padded hand grips.
- Comes with 5" wheels.
- Limited Lifetime Warranty.

**STAND ASSIST
AID WITH
RETRACTABLE
POLES!**



**LIFTWALKER WILL PROVIDE CONVENIENT SAFE
ASSISTED INDEPENDENCE ANYWHERE,
ANYTIME IT'S NEEDED!**



SPECIFICATIONS

Item #	Width (inside hand grip)	Width (inside back legs)	Overall Width	Depth (opened at base)	Height (highest setting)	Height (lowest setting)	Product Weight	Carton Shipping Weight	Weight Limit
LWWSW	19"	21.5"	25"	22.5"	39.5"	33.5"	11 lbs.	14 lbs.	300 lbs.

walkers

To order call toll free: 877.224.0946

drive™