



12026

## Locking Elevated Toilet Seat

Available  
**1 Pack**

12026

Without Arms,  
(Will become RTL12026 in 2011), 1/cs

### FEATURES:

- Designed for individuals who have difficulty sitting down or standing up from the toilet.
- Heavy-Duty molded plastic construction provides additional strength and durability.
- Improved locking device is engineered with larger, heavy-duty "worm screw" and locking plate to ensure a safe, secure fit.
- Lightweight and portable.
- Fits most toilets.
- No tools required for installation.
- Easy to clean.
- Limited Lifetime Warranty.



12014

## Premium Plastic Elevated, Regular/Elongated Toilet Seat with Lock

Available  
**1 Pack**

12014

Without Arms, 1/cs

### FEATURES:

- Designed for individuals who have difficulty sitting down or standing up from the toilet.
- Wide opening in front and back for personal hygiene.
- Lightweight and portable.
- Standard locking mechanism ensures safety.
- Fits most elongated toilets.
- No tools required for installation.
- Easy to clean.
- Limited Lifetime Warranty.



12063

12062

## Raised Toilet Seat with/without Lid

Available  
**1 Pack**

12062

Without Lid, 2", (Will become RTL12062 in 2011), 1/cs

12063

With Lid, 2", (Will become RTL12063 in 2011), 1/cs

12064

Without Lid, 4", (Will become RTL12064 in 2011), 1/cs

12065

With Lid, 4", (Will become RTL12065 in 2011), 1/cs

12066

Without Lid, 6", (Will become RTL12066 in 2011), 1/cs

12067

With Lid, 6", (Will become RTL12067 in 2011), 1/cs

12062-C

Replacement Locking Clamps, 2/bx.

### FEATURES:

- Designed for individuals who have difficulty sitting down or standing up from the toilet.
- Two hygiene cutouts provide maximum convenience.
- Designed for use with existing toilet seat.
- Lightweight and portable.
- Easily attaches to toilet bowl and locks in place with two rear locks.
- No tools required for installation.
- Easy to clean.
- Fits most toilets.
- Limited Lifetime Warranty.



**LOCKS  
SECURELY IN  
PLACE WITH  
LOCKING  
KNOB &  
CLAMP!**

**COMPARE TO:  
SAVANNAH®**

### SPECIFICATIONS

Item #	Width (seat)	Depth (seat)	Seat Height	Product Weight	Carton Shipping Weight	Weight Limit
12026	17"	16.5"	5"	4.5 lbs.	6 lbs.	300 lbs.
12014	17"	16.5"	4.5"	5 lbs.	7 lbs.	300 lbs.
12062/3	14"	16"	2"	2/3 lbs.	4/5 lbs.	350 lbs.
12064/5	14"	16"	4"	2/3 lbs.	4/5 lbs.	350 lbs.
12066/7	14"	16"	6"	3.5/4.5 lbs.	6/7.5 lbs.	350 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

bath safety

To order call toll free: 877.224.0946

drive™