



12034

**TOOLS INCLUDED FOR INSTALLATION**



**Tub Rail**

- 12034 8" x 11", 1/bx
- 12035\* 7" x 17", 1/bx
- 12043\* 15" x 17", 1/bx

**FEATURES:**

- Durable, chrome-plated steel is attractive and easy-to-clean.
- Rubber pads on all tub contact points eliminate scratches.
- Knurling around grab bar diameter provides a no-slip grip.
- Attractive retail packaging.
- Easy to install (not for use on fiberglass tubs).
- Limited Lifetime Warranty.

**Clamp on Tub Rail**

- 12036-ADJ Adjustable Height, White, 1/bx
- 12036 White, Non-Adjustable, 1/bx
- 12041 White, Non-Adjustable, 1/bx (Rail is perpendicular to tub wall)

**FEATURES:**

- Durable, white powder-coated steel construction is attractive and easy-to-clean.
- 14" - 17" Height Adjustment (for 12036-ADJ only).
- Handle design provides secure hand position.
- Width adjusts from 3" - 7".
- Stainless steel locking mechanism and non-wearing rubber pads provide scratch proof security.
- Fits most tubs (not for use with fiberglass tubs).
- Tool free installation.
- Attractive retail packaging.
- Limited Lifetime Warranty.



12036-ADJ



**14"-17" HEIGHT ADJUSTMENT**



**TOOL FREE INSTALLATION**

**SPECIFICATIONS**

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

Item #	Adjustable Clamp (smallest)	Adjustable Clamp (widest)	Height (lowest to highest)	Product Weight	Carton Shipping	Weight Limit
12034	3.5"	5.5"	not adjustable	4.4 lbs.	6.2 lbs.	300 lbs.
12035	3.5"	5.5"	not adjustable	5.1 lbs.	6.4 lbs.	300 lbs.
12043	3.75"	5.5"	not adjustable	5.6 lbs.	7.6 lbs.	300 lbs.
12036-ADJ	3"	7"	14" - 17"	6.8 lbs.	8.9 lbs.	300 lbs.
12036	3"	7"	not adjustable	6.4 lbs.	8.4 lbs.	300 lbs.
12041	3"	7"	not adjustable	6.4 lbs.	8.4 lbs.	300 lbs.

bath safety  
To order call toll free: 877.224.0946

**drive**  
WWW.DRIVEMEDICAL.COM