



COMPARE TO:
drive Brewer® Clinton®
 13071 43465 M20

Mayo-Instrument Stand, with 5 Casters

13071 1/cs
 13035T Replacement Tray, 1/cs
 (For use with 13035, 13045
 and 13071)

FEATURES:

- Casters provide additional mobility and convenience.
- Removable stainless steel tray measures 19" x 12¹/₂".
- Tray height adjusts from 29¹/₂" to 47" with lock.
- Limited Lifetime Warranty.



COMPARE TO:
drive Brewer® Clinton®
 13035 43460 M21

Mayo-Instrument Stand, Single Post

13035 1/cs
 13035T Replacement Tray, 1/cs
 (For use with 13035 and 13045)

FEATURES:

- Single post stand with "U" base design provides support and versatility.
- Removable stainless steel tray measures 19" x 12⁵/₈".
- Tray height adjusts from 32¹/₈" to 50" with lock.
- Two, 3" casters provide additional mobility.
- Limited Lifetime Warranty.



COMPARE TO:
drive Brewer® Clinton®
 13045 11100 M22

Mayo- Instrument Stand, Double Post

13045 1/cs
 13035T Replacement Tray, 1/cs
 (For use with 13035 and 13045)

FEATURES:

- Double post provides added support for use in surgical procedures.
- Removable stainless steel tray measures 19" x 12⁵/₈".
- Adjustable height from 34" to 54¹/₂".
- Two hooded casters.
- Limited Lifetime Warranty.

SPECIFICATIONS

Item #	Table/Tray Width	Table/Tray Depth	Height (highest)	Height (lowest)	Base Width	Base Depth	Product Weight	Carton Shipping Weight
13071	19"	12.5"	47"	29.5"	21"	21"	11 lbs.	14 lbs.
13035	19"	12.625"	49"	32"	22"	18.5"	9 lbs.	12 lbs.
13045	19"	12.625"	54.5"	34"	21"	19"	11 lbs.	14 lbs.

Due to product improvements and changes, all images in this catalog are subject to change.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

patient room

drive
 WWW.DRIVEMEDICAL.COM

To order call toll free: 877.224.0946