



13079

Revolving, Pneumatic Adjustable Height Stool with Plastic Base

13079 Black, 1/cs

FEATURES:

- Standard plastic base.
- Pneumatic adjustment provides easy-touch height adjustment from 14.5" to 20".
- Solid steel shaft.
- 14" diameter round seat.
- 4" thick padded seat (black).
- Black, composite base.
- 5 Casters.
- 250 lb. Weight Capacity.
- Limited Lifetime Warranty.

COMPARE TO:		
drive	Brewer®	Clinton®
13079	22500	2130



13078

Revolving, Pneumatic Adjustable Height Stool with Metal Base

13078 Black, 1/cs

FEATURES:

- Standard metal base.
- Pneumatic adjustment provides easy-touch height adjustment from 19" to 22.5".
- Solid steel shaft.
- 14" diameter round seat.
- 4" thick padded seat (black).
- 5 Casters.
- 250 lb. Weight Capacity.
- Limited Lifetime Warranty.

COMPARE TO:	
drive	Brewer®
13078	11001

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.