



Raised Rails and Custom Sizes Available!



drive™ Mason™ Multi-Ply ShearCare1500 Pressure Redistribution Foam Mattress

Item #	Product Description	Product Dimensions	Product HCPCS	UOM
1500SC-1-FB	Multi-layered/Multi-zoned Bariatric Foam Mattress	76" (L) x 36" (W) x 6" (H)	E0184	1/ea
1500SC-2-FB	Multi-layered/Multi-zoned Bariatric Foam Mattress	80" (L) x 36" (W) x 6" (H)	E0184	1/ea
1500SC-3-FB	Multi-layered/Multi-zoned Bariatric Foam Mattress	84" (L) x 36" (W) x 6" (H)	E0184	1/ea
1500SC-42-FB	42" Multi-layered/Multi-zoned Bariatric Foam Mattress	80" (L) x 42" (W) x 6" (H)	E0184	1/ea
1500SC-48-FB	48" Multi-layered/Multi-zoned Bariatric Foam Mattress	80" (L) x 48" (W) x 6" (H)	E0184	1/ea
1500SC-54-FB	54" Multi-layered/Multi-zoned Bariatric Foam Mattress	80" (L) x 54" (W) x 6" (H)	E0184	1/ea

- Multi-layered pressure redistribution foam mattress specifically designed for the high-risk bariatric patient.
- The top layer is a high-resilient, high-density die-cut foam for maximum support and comfort.
- The heel section is made of Lura-Quilt™ foam which is an aerated foam that allows air to pass through keeping the patient cool and dry.
- The bottom layer is made of LX45 foam which is a high-density, high ILD foam that provides patient support and stability.
- Designated head and foot sections prevent installation errors.
- Top Cover is a 2-way stretch nylon which reduces heat buildup and provides a low-shear surface that is anti-microbial and water-resistant.
- Concealed zipper and a Barrier Stop Over-Flap prevents liquids from contaminating mattress core.
- Meets Federal Fire Code CFR 16 part 1633.

specifications	
Cover Material	Fluid-Resistant Nylon Top Waterproof Vinyl Bottom
Mattress Layers of Foam	2
Side Rails	Yes
Dedicated Heel Section	Yes
Fire Rating	CFR 16 Part 1633 Federal Fire Code
Weight Capacity	36" 450 lbs., 42" 650 lbs., 48" 750 lbs., 54" 850 lbs.
Warranty	1 Year Limited Foam 1 Year Cover