

FULL "PAN"
SLEEP SURFACE


NEW & IMPROVED
ELECTRONICS!



BARIATRIC
750 LB
WEIGHT CAPACITY

New & Improved Full Electric Bariatric Bed

- 15302** Bariatric Bed, 48" Width, 750 lb Capacity, 1/cs
- 15302BV-1HR** Bariatric Bed with 1 pair "T" Rails, (item #15021/15021BV) (Has replaced 15302-1HR)
- 15302BV-2HR** Bariatric Bed with 2 pair "T" Rails, (item #15021/15021BV) (Has replaced 15302-2HR)
- 15302BV-PKG** Bariatric Bed Package (contains Bariatric Bed, 48" Bariatric Foam Mattress, item #15310, and 1 pair "T" Rails, item #15021/15021BV) (Has replaced 15302-PKG)
- 15302BV-PKG-2** Bariatric Bed Package (contains Bariatric Bed, 48" Bariatric Foam Mattress, item #15310, and 2 pair "T" Rails, item #15021/15021BV) (Has replaced 15302-PKG-2)



***BED PACKAGES NOW
AVAILABLE WITH
BROWN VEIN RAILS!**

FEATURES:

- Quiet, smooth operation.
- Heavy duty frame ensures strength and patient safety.
- All steel construction.
- Larger sleep surface than conventional bed.
- Quiet, smooth operation.
- Split pan design with removable bed ends is easy to set-up.
- No tools necessary for set-up.
- Hand control (included) provides for multiple bed positioning for patients.
- Emergency manual crank included.
- Low air loss/alternating pressure mattress available.
- Warranty: (Original purchaser) Lifetime on welds, 5 years on frame, 1 year on all other parts and components.

Bariatric Mattress

- 15310** 48" x 80" x 7",
Weight Capacity 750 lbs., 1/cs

FEATURES:

- Top cover: Unitek 11 Nylon.
- Bottom cover: Vinyl.
- Top layer: 2.5 lb density foam.
- Middle layer: 2.5 lb density foam.
- Base: 4 lb density foam.
- Fire retardant: Compliance with 16 CRF 1633.
- Warranty: 1 year.

SPECIFICATIONS

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

Item #	Bed Height (min)	Bed Height (max)	Overall Length	Overall Width	Product Weight	Carton Shipping Weight	Weight Limit
15302	18.25"	26.25"	88"	48"	328 lbs.	348 lbs.	750 lbs.

bariatric beds

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM