

**NEW!**



**drive**  
**Disposable Canister Kit**

**(1P)** 22330 1/ea

22330-12 12/cs

HCPSCS code: A7000 + A7002

Includes: 6' suction tube, 7" and 3" suction canister tubing, hydrophobic filter

**FEATURES:**

- 800cc disposable suction canister kit
- Raised etched cc increments
- Lid comes with auto shut off float valve

**(1P)** = 1-Pack Available

respiratory

To order call toll free: 877.224.0946

**drive**<sup>TM</sup>