



drive™ Mason™ Multi-Ply ShearCare300 Pressure Redistribution Foam Mattress

Item #	Product Description	Product Dimensions	Product HCPCS	UOM
300SC-1-FB	Single Layer/Multi-zoned Foam Mattress	76" (L) x 36" (W) x 5.5" (H)	E0184	1/ea
300SC-2-FB	Single Layer/Multi-zoned Foam Mattress	80" (L) x 36" (W) x 5.5" (H)	E0184	1/ea
300SC-3-FB	Single Layer/Multi-zoned Foam Mattress	84" (L) x 36" (W) x 5.5" (H)	E0184	1/ea
300SC-1-RR-FB	Single Layer/Multi-zoned Foam Mattress w/3" Elevated Perimeter Cut-out	76" (L) x 36" (W) x 5.5" (H)	E0184	1/ea
300SC-2-RR-FB	Single Layer/Multi-zoned Foam Mattress w/3" Elevated Perimeter Cut-out	80" (L) x 36" (W) x 5.5" (H)	E0184	1/ea
300SC-3-RR-FB	Single Layer/Multi-zoned Foam Mattress w/3" Elevated Perimeter Cut-out	84" (L) x 36" (W) x 5.5" (H)	E0184	1/ea

- Pressure Redistribution mattress designed for home care and long term care.
- This mattress is divided into 3 distinct zones, comprised of 66 channel cuts for differential pressure redistribution.
- Designated head and foot sections prevent installation errors.
- Top Cover is a 2-way stretch nylon which reduces heat buildup and provides a low-shear surface that is anti-microbial and water-resistant.
- Concealed zipper and Mason's Barrier Stop Over-Flap prevents liquids from contaminating mattress core.
- Conforms to Federal Fire Code CFR 16 part 1633.

specifications	
Cover Material	Fluid-Resistant Nylon Top Waterproof Vinyl Bottom
Mattress Layers of Foam	1
Side Rails	Available
Dedicated Heel Section	Yes
Fire Rating	CFR 16 Part 1633 Federal Fire Code
Weight Capacity	300 lbs.
Warranty	3 Year Limited Foam 1 Year Cover