



Spring-Ease™ Extra-Firm Support Innerspring Mattress

Ortho-Coil™ Super-Firm Support Innerspring Mattress

Item #	Product Dimensions	Product HCPCS	UOM
3637-1SE	76" (L) x 36" (W) x 6" (H)	EO271	1/ea
3637-2SE	80" (L) x 36" (W) x 6" (H)	EO271	1/ea
3637-3SE	84" (L) x 36" (W) x 6" (H)	EO271	1/ea

- Deluxe innerspring mattress provides extra-firm support.
- Multi-layered with two 1" high-density, firm, combustion modified fiber toppers on the top and one 1" high density, firm, combustion-modified fiber topper on the bottom.
- The Masongard vinyl waterproof cover is anti-microbial/anti-bacterial.
- The mattress is fire-retardant and conforms to CFR 16 part 1633.

specifications

Cover Material	Fluid-Proof Vinyl
Weight Capacity	300 lbs.
Warranty	1 Year

Item #	Product Dimensions	Product HCPCS	UOM
3637-1OC	76" (L) x 36" (W) x 6" (H)	EO271	1/ea
3637-2OC	80" (L) x 36" (W) x 6" (H)	EO271	1/ea
3637-3OC	84" (L) x 36" (W) x 6" (H)	EO271	1/ea

- Super-firm innerspring mattress provides ultimate support and comfort.
- The two polyester fiber pads on the top and bottom of the mattress are perfect for heavier patients.
- The Masongard vinyl waterproof cover is anti-microbial/anti-bacterial.
- The mattress is fire-retardant and conforms to CFR 16 part 1633.

specifications

Cover Material	Fluid-Proof Vinyl
Weight Capacity	350 lbs.
Warranty	1 Year