



Custom Sizes Available  
Upon Request



**drive** Mason™ The Sierra™ Gel Overlay

Item #	Product Description	Product Dimensions	Product HCPCS	UOM
3671	Gel Overlay	78" (L) x 34" (W) x 3" (H)	E0185	1/ea
3671-B	Gel Overlay	78" (L) x 34" (W) x 4" (H)	E0185	1/ea
3671-42	Gel Overlay	78" (L) x 42" (W) x 3" (H)	E0185	1/ea
3671B-42	Gel Overlay	78" (L) x 42" (W) x 4" (H)	E0185	1/ea
3671B-48	Gel Overlay	78" (L) x 48" (W) x 4" (H)	E0185	1/ea
3671B-54	Gel Overlay	78" (L) x 54" (W) x 4" (H)	E0185	1/ea

- The Sierra incorporates a series of gel bladders encapsulated between high-density polyurethane sculptured foam and a support base.
- This combination creates a therapeutically advanced yet cost effective support surface.
- Optimal environment for prevention/treatment of pressure ulcers.
- One-Piece construction requires no assembly.
- Easy to clean, low shear cover.
- Elastic corner straps to secure overlay to mattress.
- Convoluted high-density polyurethane foam topper.
- Split bladders contain an aqueous non-toxic gel solution and prevents migration.

specifications	
Cover Material	Fluid-Resistant Nylon Top   Vinyl Bottom
Weight Capacity	300 - 650 lbs.
Warranty	1 Year