



**Entry-Level, Ultra-Soft
LTC Mattress**



drive™ Mason™ Multi-Ply 6500 Lite Pressure Redistribution Foam Mattress

Item #	Product Description	Product Dimensions	Product HCPCS	UOM
6500-LT-1-FB	Multi-Layered/Multi-zoned Foam Mattress	76" (L) x 36" (W) x 6" (H)	E0184	1/ea
6500-LT-2-FB	Multi-Layered/Multi-zoned Foam Mattress	80" (L) x 36" (W) x 6" (H)	E0184	1/ea
6500-LT-3-FB	Multi-Layered/Multi-zoned Foam Mattress	84" (L) x 36" (W) x 6" (H)	E0184	1/ea
6500-LT-1-RR-FB	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter Cut-out	76" (L) x 36" (W) x 6" (H)	E0184	1/ea
6500-LT-2-RR-FB	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter Cut-out	80" (L) x 36" (W) x 6" (H)	E0184	1/ea
6500-LT-3-RR-FB	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter Cut-out	84" (L) x 36" (W) x 6" (H)	E0184	1/ea

- Multi-layered foam mattress with top layer of die-cut high-density foam in the head and torso section for excellent pressure redistribution.
- Heel section is multi-layered with independent cell technology to address this high-risk area.
- Bottom layer of high-density foam provides a firm foundation and prevents bottoming out.
- Unique SSPT Rail System™ provides solid-edge support and facilitates safe patient transfers.
- Designated head and foot sections prevent installation errors.
- Top Cover is a 2-way stretch nylon which reduces heat buildup and provides a low shear surface that is anti-microbial and water-resistant.
- Concealed zipper and a Barrier Stop Over-Flap prevents liquids from contaminating mattress core.
- Meets Federal Fire Code CFR 16 part 1633.

specifications	
Cover Material	Fluid-Resistant Nylon Top Waterproof Vinyl Bottom
Mattress Layers of Foam	2
Side Rails	Yes
Dedicated Heel Section	Yes
Fire Rating	CFR 16 Part 1633 Federal Fire Code
Weight Capacity	300 lbs.
Warranty	5 Year Limited Foam 1 Year Cover