



drive ^{Mason}™ 3" Gel/Foam Gel-U-Seat™ Cushion

Item #	Product Description	Product Dimensions	Product HCPCS	UOM
8046-1	3" Gel/Foam Cushion (Replaces Mason # 8045-1)	16" (W) x 16" (D) x 3" (H)	E2603	1/ea, 4/cs
8046-2	3" Gel/Foam Cushion (Replaces Mason # 8045-2)	18" (W) x 16" (D) x 3" (H)	E2603	1/ea, 4/cs
8046-3	3" Gel/Foam Cushion (Replaces Mason # 8045-3)	20" (W) x 16" (D) x 3" (H)	E2603	1/ea, 4/cs
8046-4	3" Gel/Foam Cushion	20" (W) x 18" (D) x 3" (H)	E2603	1/ea, 4/cs
8046-5	3" Gel/Foam Cushion	22" (W) x 18" (D) x 3" (H)	E2604	1/ea, 4/cs
8046-6	3" Gel/Foam Cushion	22" (W) x 20" (D) x 3" (H)	E2604	1/ea, 4/cs
8046-7	3" Gel/Foam Cushion	24" (W) x 18" (D) x 3" (H)	E2604	1/ea, 4/cs
8046-8	3" Gel/Foam Cushion	24" (W) x 20" (D) x 3" (H)	E2604	1/ea, 4/cs
8046-9	3" Gel/Foam Cushion	26" (W) x 18" (D) x 3" (H)	E2604	1/ea, 4/cs
8046-10	3" Gel/Foam Cushion	26" (W) x 20" (D) x 3" (H)	E2604	1/ea, 4/cs

- The Gel-U-Seat™ cushion is specifically designed to redistribute pressure to assist in the prevention of pressure ulcers while providing patient comfort.
- Cushion incorporates high-density foam encapsulating a dual-chamber gel bladder.
- Water-resistant, low-shear stretch nylon cover with tie straps to secure to wheelchair.

specifications	
Cover Material	Half P/U, Half Vinyl
Weight Capacity	Up to 20" 275 lbs., 22" and up 500 lbs.
Warranty	1 Year Limited