

Digital Dual Channel 5 Mode TENS with Two Strength Duration Modes

AGF-5X 1/bx

FEATURES:

- Digital wave technology provides a highly efficient delivery of the wave signal.
- Dual isolated channels.
- Intensity control is adjustable to 80 mA, 500 ohm load.
- Pulse frequency adjustable from 0.5-150 Hz.
- Pulse width is adjustable from 30-260 micro-seconds.
- Wave form: Asymmetrical, bi-phase, square pulse.
- Modulation selector offers 5 modes of wave delivery, Normal, Burst and 3 Modulation modes.
- Metal belt clip.
- Comes complete with carrying case, 2 lead wires, 9-volt battery, 1 package of 4 pre-gelled electrodes and operating manual.
- 3 Year Limited Warranty.



AGF-5X

Digital Electronic Muscle Stimulator with Timer

AGF-6X 1/bx

FEATURES:

- Electrical Muscle Stimulation (EMS) is a proven, efficient way of treating muscle injuries by the transmission of electronic pulses that causes passive exercise. The stimulation, contraction and relaxation cycle of EMS loosens muscle fibers, increases blood flow and stimulates muscle growth.
- Digital wave technology provides a highly efficient delivery of the wave signal.
- Dual isolated channels.
- Pulse Amplitude: 0-90mA peak to peak into 500 ohm load.
- Pulse frequency adjustable from 1-150 Hz.
- Pulse width is adjustable from 30 uS ~ 260 uS, 10 uS/step.
- Contraction time is adjustable from 1-30 seconds.
- Relaxation time is adjustable from 1-30 seconds.
- Ramp time is adjustable from 1-6 seconds.
- Metal belt clip.
- Comes complete with carrying case, 2 lead wires, 9-volt battery, 1 package of 4 pre-gelled electrodes and operating manual.
- 3 Year Limited Warranty.



AGF-6X

electrotherapy

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.