

# IV Pole

CE 1020 For use with all  
Wenzelite Safety Rollers, 1/bx

Telescoping IV pole with mounting bracket



# Oxygen Holder

CE 1250 For use with all  
Wenzelite Safety Rollers, 1/bx

Oxygen holder for a “D” or “E” cylinder

# Dial-A-Speed Tabs

CE 1400 For use with all Wenzelite Safety Rollers, 1pr/bx

Dial-A-Speed tabs provide variable resistance on one or both rear wheels. This serves as a partial brake and allows the user to control the speed of the Safety Roller. Each side is independently adjusted and, therefore, allows for alignment of asymmetrical gait patterns.



# Swivel Wheel Locking Brackets

CE 1500 For use with all Wenzelite Safety Rollers, 1pr/bx

Locking brackets for front swivel wheels ensure straight walking. They are easily unlocked for increased maneuverability. To combine the advantages of control and easy maneuverability, you can lock one side and unlock the other side.



427

mobility  
safety rollers  
accessories

Wenzelite Re/hab

To order call toll free: 877.224.0946

drive  
WWW.DRIVEMEDICAL.COM

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.