



Bed Renter II Densified Fiber Mattress

Bed Renter Foam Mattress

Item #	Product Dimensions	Product HCPCS	UOM
3500-II	76" (L) x 36" (W) x 5.5" (H)	E0272	1/ea
3501-II	80" (L) x 36" (W) x 5.5" (H)	E0272	1/ea
3502-II	84" (L) x 36" (W) x 5.5" (H)	E0272	1/ea

Item #	Product Dimensions	Product HCPCS	UOM
M3500	76" (L) x 36" (W) x 5.5" (H)	E0272	1/ea
M3501	80" (L) x 36" (W) x 5.5" (H)	E0272	1/ea
M3502	84" (L) x 36" (W) x 5.5" (H)	E0272	1/ea

- The densified fiber core is designed specifically for use in Home Care.
- Lighter weight than innerspring mattress and reduces shipping costs.
- Durable for years of practical use.
- The Masongard vinyl cover is waterproof and anti-microbial.
- Meets Federal Fire Code CFR 16 part 1633.

- The comfortable foam core is designed specifically for use in Home Care.
- Durable for years of practical use.
- The Masongard vinyl cover is waterproof and anti-microbial.
- Meets Federal Fire Code CFR 16 part 1633.
- **Replaces Mason # 3500, 3501 and 3502.**

specifications

Cover Material	Fluid-Proof Vinyl
Weight Capacity	250 lbs.
Warranty	1 Year

specifications

Cover Material	Fluid-Proof Vinyl
Weight Capacity	250 lbs.
Warranty	1 Year