



**NEW!**



**Nitro Rollator**


- 1P** RTL10266 Red, 1/cs
- 1P** RTL10266WT White, 1/cs
- 1P** RTL10266BK Black, 1/cs

**FEATURES:**

- Attractive, Euro-style design.
- Brake cable inside frame for added safety.
- Handle height easily adjusts with unique push button.
- Caster fork design enhances turning radius.
- Large 10" front casters for easy maneuverability and indoor or outdoor use.
- Easily folds with one hand to ultra-compact size for storage.
- Plastic clip holds the Nitro closed when frame is folded.
- Lightweight, aluminum frame.
- Back support height easily adjusts with tool-free thumb screw.
- Seat with fold-up handle is durable and comfortable.
- Cross brace design allows for side-to-side folding and added stability.
- Frame can be folded with one hand for ultra compact storage and portability.
- Handsome and removable zippered storage bag. Unique bag attachment keeps bag securely in place when rollator is open or folded.
- Comes in retail box.
- Limited Lifetime Warranty.



**FRAME CAN BE FOLDED WITH ONE HAND FOR ULTRA COMPACT STORAGE AND PORTABILITY!**



- A. Back support height easily adjusts with tool-free thumb screw.
- B. Handle height easily adjusts with unique push button.



**PARTS DIAGRAMS ON PAGE 32**

**1P = 1-Pack Available**

**SPECIFICATIONS**

Item #	Handle Height Highest	Handle Height Lowest	Unit Width	Unit Length	Seat (depth)	Seat (width)	Seat (height)	Product Weight	Carton Shipping Weight	Weight Limit
RTL10266 Series	38.25"	33.5"	23"	27.75"	8"	18"	20.5"	17.5 lbs.	21.5 lbs.	300 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.