



12027RA

18" WIDTH BETWEEN FRONT ARMS!

REMOVEABLE ARMS PINLOCK INTO METAL SLEEVES!

2 in 1 Locking Elevated Toilet Seat with Tool Free Removable Arms

12027RA With Removable Arms, Retail Packaged, (Has replaced 12027), 1/cs

12027RA-4BULK With Removable Arms, 4/cs (Each toilet seat individually polly bagged, 4 per carton)

FEATURES:

- Tool Free Arms can be removed or added as needed.
 - Full length padded removeable arms.
 - Arm sleeves and receivers are made of metal for more durability.
 - Heavy-duty molded plastic construction provides additional strength and durability.
 - Improved locking device is engineered with larger, heavy-duty "worm screw" and locking plate to ensure a safe, secure fit. Fits most toilets.
- Color: White.
Limited Lifetime Warranty.

REMOVEABLE ARMS ALLOW YOU TO STOCK ONE PRODUCT AND FILL THE NEEDS OF TWO.

18" WIDTH BETWEEN ARMS!



12008KDR



TOOL FREE REMOVABLE & ARMS!

K.D. 2 in 1 Locking, Elevated Toilet Seat with Tool Free Removable Arms

12008KDR (Has replaced 12008), 1/cs



FEATURES:

- Arms can be removed or added as needed.
 - Padded arms can be easily removed without tools.
 - Removable arms reduce the cost of shipping.
 - Heavy duty molded plastic construction provides additional strength and durability.
 - Clamping mechanism ensures secure locking onto toilet.
 - Prevents shifting of seat while transferring.
- Color: White
Limited Lifetime Warranty.



REMOVEABLE ARMS

SPECIFICATIONS

Item #	Width (seat)	Depth (seat)	Seat Height	Width Between Arms	Product Weight	Carton Shipping Weight	Weight Limit
12027RA	17"	16.5"	5"	18"	5.1 lbs.	8.5 lbs.	300 lbs.
12027RA-4BULK	17"	16.5"	5"	18"	5.1 lbs.	19.4 lbs.	300 lbs.
12008KDR	15.5"	16"	5"	18"	4.5 lbs.	6.5 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

bath safety

drive
WWW.DRIVEMEDICAL.COM

To order call toll free: 877.224.0946